



1,000 µl XL graduated TipOne® Filter Tip, Natural, Refills (sterile) Case



Product Highlights

- Graduated at 100 µl, 250 µl, 500 µl and 1000 µl
- Quality: Highly effective HDPE filter
- Safe: Filter contains no additives to interfere with your samples
- Retrieval: The filter will not absorb samples, allowing recovery if necessary
- Certified: Sterile, and RNase, DNase, DNA and Pyrogen free
- Sterile

Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	S1122-1730-C
Alternative Product	pf-sl-920414
Pack Size	3,840 Tips (40 Refills × 96 Tips)
Color Code Wafer	Blue
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume range	100 - 1,000 µl
Volume	1000 µl
Tip Type	XL
Material	Polypropylene
Graduations	Yes
Color	Natural
For use with	Abimed : BioPette-E, Discovery Comfort, LABMATE Optima Biohit : ePet, mLine, Proline BRAND : Transferpette, Transferpette Electronic, Transferpette S Capp : CappAero Comfort, CappAero Universal, Cappronic Single, ecopipette, Singlechannel Eppendorf : Reference, Reference 2, Research, Research plus, Research pro, Xplorer Gilson : Pipetman Neo, Pipetman P, Pipetman Ultra Hamilton : SoftGrip Pipette, SofTouch HTL : Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette Nichiryo : 31N, 50F, 50V, NLE, NPLO, NPX Rainin : Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus Sartorius : ePet, mLine, Proline, Proline Plus, Tacta Socorex : Acura 815 manual, Acura 825, Acura electro 925, Calibra digital STARLAB : ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E Thermo Fisher Scientific Inc. : Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 1,000 µl XL graduated TipOne® Filter Tip, Natural, Refills (sterile) Case

Nothing goes through!

Protect your samples against unwanted contamination by using TipOne® filter tips. The TipOne® filter is one of the best on the market and protects against aerosol cross-contamination without trapping your valuable samples. The filters contain no additives to interfere with samples.









During pipetting it is not possible to see with the eye if fine aerosols are going into your pipette. If they are, those aerosols can contaminate other samples at a later date.

TipOne® filter tips are proven to block >99% of aerosols and are among the best filters on the market. The results, shown right, speak for themselves!

**The comparative tests were conducted in an independent test institute with filter tips from other leading manufacturers.*



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1,000 µl XL graduated TipOne® Filter Tip, Natural, Refills (sterile) Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1122-1730
 1,000 µl XL graduated TipOne® Filter Tip, Natural, Refills (sterile) Case Pack Size: 3,840 Tips (40 Refills × 96 Tips)  Create PDF Data Sheet	3,840 Tips (40 Refills × 96 Tips)	S1122-1730-C
 1,000 µl XL graduated TipOne® Filter Tip, Natural, Racks (sterile) Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1122-1830
 1,000 µl XL graduated TipOne® Filter Tip, Natural, Racks (sterile) Case Pack Size: 3,840 Tips (40 Racks × 96 Tips)  Create PDF Data Sheet	3,840 Tips (40 Racks × 96 Tips)	S1122-1830-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk