



## Refill Tips for Beckman Coulter® (Biomek® FX, Biomek® NX, Biomek®3000, Multimek®)



### Product Highlights

- Tips suitable for Biomek® i-Series • Biomek® FX • Biomek® NX • Biomek® 3000 • Biomek® 4000 • Multimek®
- Space saving refill stacks that take up 40% less space
- Certified RNase, Dnase, DNA & PCR Inhibitor free
- No special set up required! Just rack up and they're ready to go

## Product Applications

PCR Setup

Assay-Setup

Serial Dilutions

Distributing Reagents

Sample or Reagent Transfer

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	See variations
Filter tip	Non-filter
For use with (SL)	Beckman Coulter®

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Refill Tips for Beckman Coulter® (Biomek® FX, Biomek® NX, Biomek®3000, Multimek®)

Robotic tip refills, suitable for Biomek® i-Series • Biomek® FX • Biomek® NX • Biomek® 3000 • Biomek® 4000 • Multimek®.  
Refills take up 40% less space than racked tips and save on plastic waste

Engineered to the same specifications as the original equipment manufacturers' tips to ensure compatibility and a consistent fit. Tips are tested on the original robotic equipment so there is no doubt they will perform to the same requirements as the disposables from the robot manufacturer themselves.









Certified free of RNase, DNase, DNA and PCR Inhibitors. Tested Pyrogen free.

Quality control is the essential factor which gives you the freedom to walk away from your automated workstation and get on with other tasks. STARLAB robotic tips are put through stringent and thorough testing. Every lot of tips is subject to an array of demanding tests, from the geometry of the tip itself to automation-specific requirements such as clean ejection from the liquid handling head. Every lot is required to meet or exceed the parameters necessary to give you that confidence to walk away.

No special set-up required! Starlab robotic tips require no instrument hardware or software adjustment, thereby providing a simple, safe and direct high-performance substitute.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">20µl Tip for Beckman Coulter®, Refill</a> <b>Volume:</b> 20 µl <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	E1076-4700
 <a href="#">20µl Tip for Beckman Coulter®, Refill (Case Size)</a> <b>Volume:</b> 20 µl <b>Sterile:</b> No <b>Pack Size:</b> 4,800 Tips (50 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	4,800 Tips (50 Refills × 96 Tips)	E1076-4700-C
 <a href="#">250µl Tip for Beckman Coulter®, Refill</a> <b>Volume:</b> 250 µl <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	E1076-0700
 <a href="#">250µl Tip for Beckman Coulter®, Refill (Case Size)</a> <b>Volume:</b> 250 µl <b>Sterile:</b> No <b>Pack Size:</b> 4,800 Tips (50 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	4,800 Tips (50 Refills × 96 Tips)	E1076-0700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)