



Robotic Reservoirs



Product Highlights

- All robotic reservoirs are certified RNase, DNase and Pyrogen free

Product Applications

PCR Setup

Assay-Setup

Serial Dilutions

Distributing Reagents

Sample or Reagent Transfer

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|-----------|---------------------------------------|
| Art. No. | See variations |
| Sterile | No |
| Pack Size | 25 Reservoirs (1 Box × 25 Reservoirs) |
| Material | Polypropylene |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Robotic Reservoirs

High quality for high throughput!





Starlab reservoirs conform to the standard SBS footprint for compatibility with most robotic systems. They have virtually no dead volume and are available in high- or lowprofile and two well styles. The sturdy polypropylene construction ensures rigidity during automation and provides excellent chemical resistance.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|---------------------------------------|------------|
|  240 ml Reservoir, Single-Well 96 Pyramid Bottoms High Profile Bottoms: Pyramid Number of troughs: 1  Create PDF Data Sheet | 25 Reservoirs (1 Box × 25 Reservoirs) | E2999-2496 |
|  84 ml Reservoir, Multi-Well 12-Channel High Profile Bottoms: Trough-shape Number of troughs: 12  Create PDF Data Sheet | 25 Reservoirs (1 Box × 25 Reservoirs) | E2999-8412 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk