



ErgoOne® Single-Channel Starter Packs



Product Highlights

- Three or four ErgoOne® single-channel pipettes
- Three TipOne® racks with compatible tips
- Three or four Self-Adhesive Shelf Holder for one Pipette
- One STARLAB Lab Timer
- Please select the volume range for more detailed info about the content of each Starter Pack.

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

| | |
|----------------|----------------|
| Art. No. | See variations |
| Volume range | 2 - 1,000 µl |
| Volume | 1000 µl |
| Channel type | 1-channel |
| Operating mode | Manual |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about ErgoOne® Single-Channel Starter Packs







Every day, pipettes are confronted with multiple tasks that must be exact and precise. Reliability is key in the laboratory and your pipette must be trusted to perform well, day in, day out. Start your new pipetting experience with ErgoOne® single-channel pipettes and original TipOne® tips - the perfect pipetting partnership! Starter Packs also include one Self-Adhesive Shelf Holder for each pipette, plus a STARLAB Lab Timer.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------|-----------------|
|  ErgoOne® Starter Pack 4 (2.5) Volume range: 0.1 - 1,000 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | STARTERPACK42.5 |
|  ErgoOne® Starter Pack 4 (10) Volume range: 0.5 - 1,000 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | STARTERPACK410 |
|  ErgoOne® Starter Pack Volume range: 20 - 1,000 µl  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | STARTERPACK |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk