



## ErgoOne® Single-Channel Starter Packs



### Product Highlights

- Two ErgoOne® single-channel pipettes: 2.5µl & 10µl
- Two racks of 10/20ul XL Graduated TipOne® Tips, (Sterile)
- Two Self-Adhesive Acrylic Pipette Shelf Holders
- One STARLAB Lab Timer

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	See variations
Channel type	1-channel
Volume range	2 - 1,000 µl
Volume	1000 µl

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about ErgoOne® Single-Channel Starter Packs







Every day, pipettes are confronted with multiple tasks that must be exact and precise. Reliability is key in the laboratory and your pipette must be trusted to perform well, day in, day out. Start your new pipetting experience with ErgoOne® single-channel pipettes and original TipOne® tips - the perfect pipetting partnership! Starter Packs also include one Self-Adhesive Shelf Holder for each pipette, plus a STARLAB Lab Timer.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">ErgoOne® Starter Pack 4 (2.5)</a> Volume range: 0.1 - 1,000 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	STARTERPACK42.5
 <a href="#">ErgoOne® Starter Pack 4 (10)</a> Volume range: 0.5 - 1,000 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	STARTERPACK410
 <a href="#">ErgoOne® Starter Pack</a> Volume range: 20 - 1,000 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	STARTERPACK

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)