



## ErgoOne® E Multi-Channel Pipettes



### Product Highlights

- Easy to use, intuitive menu structure
- Lightweight, ergonomic design
- Adjustable finger rest for individual comfort
- Significantly reduced tip loading and ejection forces
- Universal tip fit
- Reliable
- Autoclavable pipette shaft
- 12-month warranty

## Product Applications

Pipetting

Reverse Pipetting of "tricky" liquids (e. g. with high viscosity)

Mixing

Dispense

Electrophoresis (excl. 1250 µl models)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Autoclavable	Autoclavable pipette shaft
Operating mode	Electronic
Pipetting type	Air-cushion
Operation	Battery operated piston stroke / manual tip ejection

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about ErgoOne® E Multi-Channel Pipettes













### Reliable and ergonomic!

[ErgoOne](#)® E electronic [pipettes](#) provide reliable pipetting every day. Ergonomic and lightweight, [ErgoOne](#)® E [pipettes](#) are perfect for comfortable pipetting during long periods of use. [ErgoOne](#)® E and Starlab's [TipOne](#)® [pipette tips](#) provide the perfect pipetting partnership, but the optimized design of the tip cone provides the user with a convenient universal tip fit.

Four operating modes provide standard and reverse pipetting, as well as mixing and dispensing. A unique “Electrophoresis” mode (available on all models excluding 1000 µl, 1250 µl, and 5000 µl,) is also available to help you load gels easily. The long-lasting battery can also be charged during use so you are never without your pipette.











## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">ErgoOne® E 8-Channel Electronic Pipette, 0.5 – 10µl</a> <b>Channel type:</b> 8-channel <b>Volume range:</b> 0.5 - 10 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9008-0010
 <a href="#">ErgoOne® E 8-Channel Electronic Pipette, 1.0 – 20µl</a> <b>Channel type:</b> 8-channel <b>Volume range:</b> 1.0 - 20 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9008-0020
 <a href="#">ErgoOne® E 8-Channel Electronic Pipette, 5.0 – 100µl</a> <b>Channel type:</b> 8-channel <b>Volume range:</b> 5.0 - 100 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9008-0100
 <a href="#">ErgoOne® E 8-Channel Electronic Pipette, 15 – 300µl</a> <b>Channel type:</b> 8-channel <b>Volume range:</b> 15 - 300 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9008-0300
 <a href="#">ErgoOne® E 8-Channel Electronic Pipette, 50 – 1250µl</a> <b>Channel type:</b> 8-channel <b>Volume range:</b> 50 - 1,250 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9008-1250
 <a href="#">ErgoOne® E 12-Channel Pipette Electronic, 5.0 – 100µl</a> <b>Channel type:</b> 12-channel <b>Volume range:</b> 5.0 - 100 µl  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9012-0100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)











PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">ErgoOne® E 12-Channel Electronic Pipette, 0.5 – 10µl</a> <b>Channel type:</b> 12-channel <b>Volume range:</b> 0.5 - 10 µl   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9012-0010
 <a href="#">ErgoOne® E 12-Channel Electronic Pipette, 1.0 – 20µl</a> <b>Channel type:</b> 12-channel <b>Volume range:</b> 1.0 - 20 µl   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9012-0020
 <a href="#">ErgoOne® E 12-Channel Electronic Pipette, 15 – 300µl</a> <b>Channel type:</b> 12-channel <b>Volume range:</b> 15 - 300 µl   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9012-0300
 <a href="#">ErgoOne® E 12-Channel Electronic Pipette, 50 – 1250µl</a> <b>Channel type:</b> 12-channel <b>Volume range:</b> 50 - 1,250 µl   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	G9012-1250

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Replacement Battery for ErgoOne® E</a>	1 Piece (1 Box × 1 Piece)	G9999-1005
 <a href="#">Silicone Oil for Piston (for 5000 µl Single-Channel and all Multi-Channel Pipettes)</a>	1 Piece (1 Pack × 1 Piece)	G9999-1007
 <a href="#">Fluorstatic Grease for Guiding Axis (for Multi-Channel Pipettes)</a>	1 Piece (1 Pack × 1 Piece)	G9999-1008
 <a href="#">Individual Stand (for all ErgoOne® E Multi-Channel Pipettes up to 300 µl)</a>	1 Piece (1 Box × 1 Piece)	G9999-1004
 <a href="#">Shelf Holder for one ErgoOne® E pipette</a>	1 Piece (1 Box × 1 Piece)	G9999-1009
 <a href="#">Wall Mount for one ErgoOne® E pipette</a>	1 Piece (1 Box × 1 Piece)	G9999-1010
 <a href="#">Universal stand for one ErgoOne® E pipette - all volumes</a>	1 Piece (1 Box × 1 Piece)	G9999-1011
 <a href="#">AC Adapter for Europe/UK (230V/50 Hz)</a>	1 Piece (1 Box × 1 Piece)	G9999-1000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)