



## Aspiration Pipettes (sterile), Ind. wrapped



### Product Highlights

- Non-cytotoxic and non-haemolytic USP Class VI polystyrene
- Gamma sterilized to SAL  $10^{-6}$  (ISO 11137)
- Certified RNase, DNase, latex and pyrogen free
- Individually wrapped in traceable paper-plastic peel-back packs showing lot number and expiry date
- Specially designed carton with perforated tear-strip doubles as efficient, bench top dispenser
- ISO 9001 and ISO 13485 manufactured

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	See variations
Sterile	Yes
Material	Polystyrene
For use with (SL)	Aspirator

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about **Aspiration Pipettes (sterile), Ind. wrapped**






Sterile, polystyrene pipette without a filter or graduations. Aspiration pipettes are used for vacuum aspirating procedures or transfers of fluids where accuracy is not a requirement.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">2.0 ml Aspiration Pipette (sterile), Ind. wrapped</a> <b>Volume:</b> 2.0 ml  <a href="#">Create PDF Data Sheet</a>	500 Individually wrapped sachets (1 Box × 500 Individually wrapped sachets)	E4861-0002
 <a href="#">5.0 ml Aspiration Pipettes (sterile), Ind. wrapped</a> <b>Volume:</b> 5.0 ml  <a href="#">Create PDF Data Sheet</a>	200 Individually wrapped sachets (1 Box × 200 Individually wrapped sachets)	E4861-0005
 <a href="#">10.0 ml Aspiration Pipettes (sterile), Ind. wrapped</a> <b>Volume:</b> 10.0 ml  <a href="#">Create PDF Data Sheet</a>	200 Individually wrapped sachets (1 Box × 200 Individually wrapped sachets)	E4861-0010

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)