



StarStore 81 Storage Box



Product Highlights

- 81-place polypropylene storage box with fixed 9x9 dividers for tubes up to 13mm
- Transparent lift-off lid with frosted writing panel on side
- Vent holes in base to aid draining and lid removal
- Printed reference grid on base
- Temperature use -90 °C to 121 °C
- Dimensions (WxDxH): 130 x 130 x 45 mm
- Max. tube height: 60 mm
- Autoclavable

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Dimensions (W × D × H)	130 mm × 130 mm × 45 mm
Dimensions (W × D)	130 mm × 130 mm
Max. tube length	60 mm
Autoclavable	Yes
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Material	Polypropylene
Color	Various
Places	81

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about StarStore 81 Storage Box

















81-place polypropylene storage box with fixed, 9x9 dividers for tubes up to 13mm, eg. 1.5/2.0 ml [microcentrifuge tubes](#) or screw cap microtubes. The non-hinged, lift-off lid enables tubes of varying heights to be stored inside. Standard 130x130mm footprint for use in [freezer racking](#).

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 StarStore 81, Natural Color: Natural  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5040
 StarStore 81, Neon Blue Color: Neon blue  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5041
 StarStore 81, Neon Green Color: Neon green  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5042
 StarStore 81, Neon Orange Color: Neon orange  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5043
 StarStore 81, Neon Pink Color: Neon pink  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5044
 StarStore 81, Neon Yellow Color: Neon yellow  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5046
 StarStore 81, Neon Mixed Color: Neon mixed  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5048
 StarStore 81, Black Color: Black  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2381-5049

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk