



## StarStore 100 Storage Box



### Product Highlights

- 100-place polypropylene storage box with fixed 10x10 dividers for tubes up to 13mm
- Hinged, lockable lid with positioning lugs for stable stacking
- Moulded reference grid on base.
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 140 x 155 x 55 mm

## Product Applications

Sample storage

Storage of microcentrifuge tubes or screw cap microtubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|                        |                          |
|------------------------|--------------------------|
| Art. No.               | See variations           |
| Pack Size              | 5 Pcs. (1 Pack × 5 Pcs.) |
| Dimensions (W × D × H) | 140 mm × 155 mm × 55 mm  |
| Dimensions (W × D)     | 140 mm × 155 mm          |
| External width         | 140 mm                   |
| Autoclavable           | Yes                      |
| Material               | Polypropylene            |
| Places                 | 100                      |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about StarStore 100 Storage Box















100-place polypropylene storage box with fixed, 10x10 dividers for tubes up to 13mm, eg. 1.5 / 2.0 ml microcentrifuge tubes or screw cap microtubes. Hinged, lockable lid with positioning lugs for stable stacking. Can be used in the freezer racking made specifically for use with this storage box.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

| PRODUCT NAME  |  | PACKAGING SIZE           | ART. NO.   |
|---|--|--------------------------|------------|
|    | <a href="#">StarStore 100-Place Storage Box, Natural</a><br><b>Color:</b> Natural<br> <a href="#">Create PDF Data Sheet</a>           | 5 Pcs. (1 Pack × 5 Pcs.) | I2310-5800 |
|    | <a href="#">StarStore 100-Place Storage Box, Neon Blue</a><br><b>Color:</b> Neon blue<br> <a href="#">Create PDF Data Sheet</a>       | 5 Pcs. (1 Pack × 5 Pcs.) | I2310-5841 |
|    | <a href="#">StarStore 100-Place Storage Box, Neon Green</a><br><b>Color:</b> Neon green<br> <a href="#">Create PDF Data Sheet</a>     | 5 Pcs. (1 Pack × 5 Pcs.) | I2310-5842 |
|   | <a href="#">StarStore 100-Place Storage Box, Neon Orange</a><br><b>Color:</b> Neon orange<br> <a href="#">Create PDF Data Sheet</a> | 5 Pcs. (1 Pack × 5 Pcs.) | I2310-5843 |
|  | <a href="#">StarStore 100-Place Storage Box, Neon Pink</a><br><b>Color:</b> Neon pink<br> <a href="#">Create PDF Data Sheet</a>     | 5 Pcs. (1 Pack × 5 Pcs.) | I2310-5844 |
|  | <a href="#">StarStore 100-Place Storage Box, Neon Yellow</a><br><b>Color:</b> Neon yellow<br> <a href="#">Create PDF Data Sheet</a> | 5 Pcs. (1 Pack × 5 Pcs.) | I2310-5846 |
|  | <a href="#">StarStore 100-Place Storage Box, Neon Mixed*</a><br><b>Color:</b> Neon mixed<br> <a href="#">Create PDF Data Sheet</a>  | 5 Pcs. (1 Pack × 5 Pcs.) | I2310-5848 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)