



100-Place Long Storage Box



Product Highlights

- 100-place polypropylene storage box with fixed 5x20 dividers for tubes up to 13mm, eg 1.5/2.0ml microcentrifuge tubes or screw cap microtubes
- Lift-off lid with moulded matrix
- Temperature use: -90°C to 121°C
- Autoclavable
- Dimensions: 76 x 287 x 36 mm
- Max. tube height: 50 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Dimensions (W × D × H)	76 mm × 287 mm × 36 mm
Dimensions (W × D)	76 mm × 287 mm
Max. tube length	50 mm
External width	76 mm
Autoclavable	Yes
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Material	Polypropylene
Color	Various
Operating temperature	-90°C to 121°C
Places	100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 100-Place Long Storage Box

100-place polypropylene storage box with fixed, 5x20 dividers for tubes up to 13mm. Lift-off lid with moulded matrix.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 100-Place Long Storage Boxes, Natural Color: Natural  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	E2100-1100
 100-Place Long Storage Boxes, Blue Color: Blue  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	E2100-1101-A
 100-Place Long Storage Boxes, Green Color: Green  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	E2100-1102
 100-Place Long Storage Boxes, Red Color: Red  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	E2100-1104

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk