



## 100-Place Long Storage Box



### Product Highlights

- 100-place polypropylene storage box with fixed 5x20 dividers for tubes up to 13mm, eg 1.5/2.0ml microcentrifuge tubes or screw cap microtubes
- Lift-off lid with moulded matrix
- Temperature use: -90°C to 121°C
- Autoclavable
- Dimensions: 76 x 287 x 36 mm
- Max. tube height: 50 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	See variations
Material	Polypropylene
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Operating temperature	-90°C to 121°C
Places	100
Autoclavable	Yes
Max. tube length	50 mm
External width	76 mm
Dimensions (W × D × H)	76 mm × 287 mm × 36 mm
Dimensions (W × D)	76 mm × 287 mm
Color	Various

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 100-Place Long Storage Box









100-place polypropylene storage box with fixed, 5x20 dividers for tubes up to 13mm. Lift-off lid with moulded matrix.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">100-Place Long Storage Boxes, Natural</a> <b>Color:</b> Natural  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	E2100-1100
 <a href="#">100-Place Long Storage Boxes, Blue</a> <b>Color:</b> Blue  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	E2100-1101-A
 <a href="#">100-Place Long Storage Boxes, Green</a> <b>Color:</b> Green  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	E2100-1102
 <a href="#">100-Place Long Storage Boxes, Red</a> <b>Color:</b> Red  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	E2100-1104

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)