



## StarStore 50 (hinged)



### Product Highlights

- 50-place polypropylene storage box with fixed 10x10 dividers for tubes up to 13mm
- Hinged, lockable lid with positioning lugs for stable stacking
- Moulded reference grid on base.
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 140 x 80 x 55 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	140 mm × 80 mm × 55 mm
Dimensions (W × D)	140 mm × 80 mm
Max. tube length	13 mm
Autoclavable	Yes
Pack Size	5 Pcs. (5 Pcs.)
Material	Polypropylene
Color	Various
Operating temperature	-90°C to 121°C
Places	50
For use with (SL)	for tubes up to 13mm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about StarStore 50 (hinged)











50-place polypropylene storage box with fixed, 5x10 dividers for tubes up to 13mm, eg. 1.5/2.0ml [microcentrifuge](#) tubes or screw cap microtubes. Hinged, lockable lid with positioning lugs for stable stacking.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">StarStore 50 Neon (Hinged), Natural</a> Color: Natural  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5800
 <a href="#">StarStore 50 Neon (Hinged), Neon Blue</a> Color: Neon blue  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5841
 <a href="#">StarStore 50 Neon (Hinged), Neon Green</a> Color: Neon green  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5842
 <a href="#">StarStore 50 Neon (Hinged), Neon Yellow</a> Color: Neon yellow  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5846
 <a href="#">StarStore 50 Neon (Hinged), Neon Mixed</a> Color: Neon mixed  <a href="#">Create PDF Data Sheet</a>	5 Pcs. (1 Pack × 5 Pcs.)	I2350-5848

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)