



## Polycarbonate Cryogenic Storage Boxes



### Product Highlights

- Polycarbonate storage boxes for ultra-low temperature use
- 81-place has grid dividers
- 100-place has pronged dividers
- White base with transparent, lift-off lid with printed reference grid
- Drain holes in base for venting and safety
- Temperature use -196 °C to 121 °C
- Autoclavable
- Dimensions (WxDxH): 133 x 133 x 53 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	133 mm × 133 mm × 53 mm
Dimensions (W × D)	133 mm × 133 mm
External width	133 mm
Autoclavable	Yes
Material	Polycarbonate
Color	White

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Polycarbonate Cryogenic Storage Boxes





Two polycarbonate storage boxes for ultra-low temperature use. Standard 130x130 mm footprint for use in standard freezer racking.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">81-place Polycarbonate Storage Box</a> <b>Places: 81</b>	5 Pcs. (1 Pack × 5 Pcs.)	I2381-2305
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">100-place Polycarbonate Storage Box</a> <b>Places: 100</b>	5 Pcs. (1 Pack × 5 Pcs.)	I2510-5005
	 <a href="#">Create PDF Data Sheet</a>		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)