



## Centrifuge Tube Storage and Transport Box



### Product Highlights

- Sturdy polypropylene storage / transport box for centrifuge tubes
- One holds 36 x 15ml tubes; the other holds 16 x 50ml tubes
- Transparent, lift-off lid
- Base can be removed for use in waterbath
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 145 x 145 x 127 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|                        |                          |
|------------------------|--------------------------|
| Art. No.               | See variations           |
| Rack Type              | Tube Rack                |
| Material               | Polypropylene            |
| Pack Size              | 2 Pcs. (2 Pcs.)          |
| Autoclavable           | Yes                      |
| Dimensions (W × D × H) | 145 mm × 145 mm × 127 mm |
| Dimensions (W × D)     | 145 mm × 145 mm          |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Centrifuge Tube Storage and Transport Box





Two polypropylene storage / transport boxes for centrifuge tubes. One holds 36 x 15ml tubes; the other holds 16 x 50ml tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

| PRODUCT NAME   | PACKAGING SIZE          | ART. NO.   |
|--|-------------------------|------------|
| <br><a href="#">Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue</a><br><b>Places: 16</b><br> <a href="#">Create PDF Data Sheet</a>  | 2 Pcs. (1 Box × 2 Pcs.) | E9025-1650 |
| <br><a href="#">Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green</a><br><b>Places: 36</b><br> <a href="#">Create PDF Data Sheet</a> | 2 Pcs. (1 Box × 2 Pcs.) | E9025-3615 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)