



Centrifuge Tube Storage and Transport Box



Product Highlights

- Sturdy polypropylene storage / transport box for centrifuge tubes
- One holds 36 x 15ml tubes; the other holds 16 x 50ml tubes
- Transparent, lift-off lid
- Base can be removed for use in waterbath
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 145 x 145 x 127 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

| | |
|------------------------|--------------------------|
| Art. No. | See variations |
| Dimensions (W × D × H) | 145 mm × 145 mm × 127 mm |
| Dimensions (W × D) | 145 mm × 145 mm |
| Autoclavable | Yes |
| Pack Size | 2 Pcs. (2 Pcs.) |
| Rack Type | Tube Rack |
| Material | Polypropylene |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Centrifuge Tube Storage and Transport Box

Two polypropylene storage / transport boxes for centrifuge tubes. One holds 36 x 15ml tubes; the other holds 16 x 50ml tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|-------------------------|------------|
|  Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue Places: 16  Create PDF Data Sheet | 2 Pcs. (1 Box × 2 Pcs.) | E9025-1650 |
|  Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green Places: 36  Create PDF Data Sheet | 2 Pcs. (1 Box × 2 Pcs.) | E9025-3615 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk