



StarRack Multi



Product Highlights

- Impact-resistant acetal polymer rack for cryogenic storage vials.
- Excellent chemical resistance
- Autoclavable

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Dimensions (W × D × H)	293 mm × 115 mm × 65 mm
Dimensions (W × D)	293 mm × 115 mm
Autoclavable	Yes
Rack Type	Multi-Tube Rack
Material	Acetal plastic
Color	Various

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about StarRack Multi

18, 60, or 84-place acetal polymer tube racks. Convenient carrying handles on either end of rack. Ideal for use in the waterbath as the rack will sink to the bottom without toppling over! Dimensions (WxDxH): 293 x 115 x 65 mm.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 84-Place StarRack Multi for 13mm Tubes, Mixed Places: 84 Color: Mixed  Create PDF Data Sheet	5 Piece (1 Rack × 5 Piece)	E2396-9759
 18-Place StarRack Multi for 50ml Centrifuge Tubes, Mixed Places: 18 Color: Mixed  Create PDF Data Sheet	5 Piece (1 Rack × 5 Piece)	E2396-9779
 60-Place StarRack Multi for 16mm Centrifuge Tubes, Mixed Places: 60 Color: Mixed  Create PDF Data Sheet	5 Piece (1 Rack × 5 Piece)	E2396-9769

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk