



## StarRack 80



### Product Highlights

- 80-place polypropylene fraction collection rack with a 5x16 array for tubes up to 11mm, eg. 1.5/2.0ml microcentrifuge tubes
- Alpha numeric reference grid
- The rack does not have a lid, but an optional, natural polypropylene storage box with a hinged lid into which the rack fits can be used for extra security during freezer storage
- Temperature use: -90 °C to 121 °C
- Autoclavable
- Dimensions (WxDxH): 226 x 67 x 28 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	226 mm × 67 mm × 28 mm
Dimensions (W × D)	226 mm × 67 mm
External width	226 mm
Autoclavable	Yes
Rack Type	Multi-Tube Rack
Material	Polypropylene
Color	Various
Places	80

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about StarRack 80







One of the most popular racks in the lab. Polypropylene provides high chemical resistance. Temperature use: -90 °C to 121 °C. Autoclavable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">StarRack 80 polypropylene fraction collection rack, Natural</a> <b>Color:</b> Natural	5 Pcs. (1 Pack × 5 Pcs.)	E2380-5000
	<a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarRack 80 polypropylene fraction collection rack, Mixed</a> <b>Color:</b> Mixed	5 Pcs. (1 Pack × 5 Pcs.)	E2380-5008
	<a href="#">Create PDF Data Sheet</a>		
	<a href="#">StarRack 80 polypropylene fraction collection rack, Neon Mixed</a> <b>Color:</b> Neon mixed	5 Pcs. (1 Pack × 5 Pcs.)	E2380-5048
	<a href="#">Create PDF Data Sheet</a>		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)