



CryoSecure®



Product Highlights

- White rack maintains 0 °C for a minimum of 5 hours.
- Blue rack maintains -10 °C to -20 °C for up to 3 hours.
- Dimensions (WxDxH): 172 x 110 x 105 mm with lid.

Product Applications

- Sample storage
- Sample transportation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Dimensions (W × D × H)	172 mm × 110 mm × 105 mm
Dimensions (W × D)	172 mm × 110 mm
Autoclavable	No
Color	Various

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about CryoSecure®

Double-sided racks filled with ice-gel, for use with 20 x 0.5 ml tubes on one side, or 20 x 1.5 / 2.0 ml tubes on the other.





NOTE: Not autoclavable.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 CryoSecure® (White, 0 °C) Color: White  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	I2312-2720
 CryoSecure® (Blue, -20 °C) Color: Blue  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	I2312-2721

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk