




Chest Freezer Racks for Standard Storage Boxes



Product Highlights

- Stainless steel racking with spring lock to hold the boxes in place
- For use with all standard storage boxes up to 136 x 136 x 55 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Dimensions (W × D)	140 mm × 140 mm
Dimensions (W × H)	140 mm ×
Rack Type	Chest Freezer Rack

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Chest Freezer Racks for Standard Storage Boxes

These racks for chest freezers are designed for use standard storage boxes up to 136 x 136 x 55mm.

Suitable for use with STARLAB's:

- > StarStore 130 x 130 mm range of storage boxes
- > 50 mm high cardboard storage boxes
- > Cryogenic storage boxes (excl. 81-Tall)
- > StarRack 64

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 9-Place Stainless Steel Rack for Chest Freezer Places: 9  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	I2602-9000
 13-High Aluminium Rack Places: 13  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	I2612-1300

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk