




## Chest Freezer Racks for Standard Storage Boxes



### Product Highlights

- Stainless steel racking with spring lock to hold the boxes in place
- For use with all standard storage boxes up to 136 x 136 x 55 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|                    |                    |
|--------------------|--------------------|
| Art. No.           | See variations     |
| Dimensions (W × D) | 140 mm × 140 mm    |
| Dimensions (W × H) | 140 mm ×           |
| External width     | 140 mm             |
| Rack Type          | Chest Freezer Rack |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Chest Freezer Racks for Standard Storage Boxes





These racks for chest freezers are designed for use standard storage boxes up to 136 x 136 x 55mm.

Suitable for use with STARLAB's:

- › StarStore 130 x 130 mm range of storage boxes
- › 50 mm high cardboard storage boxes
- › Cryogenic storage boxes (excl. 81-Tall)
- › StarRack 64



## All Variations

| PRODUCT NAME   | PACKAGING SIZE             | ART. NO.   |
|--|----------------------------|------------|
| <br><a href="#">9-Place Stainless Steel Rack for Chest Freezer</a><br><b>Places: 9</b><br> <a href="#">Create PDF Data Sheet</a> | 1 Piece (1 Rack × 1 Piece) | I2602-9000 |
| <br><a href="#">13-High Aluminium Rack</a><br><b>Places: 13</b><br> <a href="#">Create PDF Data Sheet</a>                        | 1 Piece (1 Rack × 1 Piece) | I2612-1300 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)