



## Upright freezer racks for 15 ml and 50 ml centrifuge tubes



### Product Highlights

- Direct storage of 15 ml or 50 ml centrifuge tubes
- Saves the time of transferring the tubes to / from storage boxes.

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Rack Type	Upright Freezer Rack
Material	Stainless Steel

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)













## More informations about Upright freezer racks for 15 ml and 50 ml centrifuge tubes

Stainless Steel (SS) racks for upright freezers, for the direct storage of 15 ml or 50 ml [centrifuge tubes](#). Saves the time of transferring the tubes to / from [storage boxes](#).

Tubes not included. Image for illustrative purposes only in order to show the storage of both 15ml and 50ml tubes.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">60-Place SS Rack for 15ml Tubes</a> <b>For use with (SL):</b> 15ml centrifuge tubes <b>Places:</b> 60   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	I2615-1000
 <a href="#">160-Place SS Rack for 15ml Tubes</a> <b>For use with (SL):</b> 15ml centrifuge tubes <b>Places:</b> 160   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	I2615-2600
 <a href="#">120-Place SS Rack for 15ml Tubes</a> <b>For use with (SL):</b> 15ml centrifuge tubes <b>Places:</b> 120   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	I2615-3600
 <a href="#">78-Place Rack for 50ml Tubes</a> <b>For use with (SL):</b> 50ml centrifuge tubes <b>Places:</b> 78   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	I2650-2600
 <a href="#">60-Place SS Rack for 50ml Tubes</a> <b>For use with (SL):</b> 50ml centrifuge tubes <b>Places:</b> 60   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	I2650-3600

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)