



Cooler Containers



Product Highlights

- The snug-fitting lid minimises thawing and evaporation
- Carrying handles
- Handy pour spout for easy emptying
- Tapered sides allow stacking when not in use
- Temperature use down to -80 °C

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	See variations
Pack Size	1 Piece (1 Piece)
Autoclavable	No
Material	Polyurethane
Color	Various

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk









More informations about Cooler Containers

Containers for temperature sensitive samples

Rigid moulded polyurethane containers for chilling your temperature sensitive samples. The Cooler Containers provide excellent insulation properties for very low temperature use with ice and ice-salt slurries.



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	4 Litre Oblong Cooler Container, Blue (380 x 250 x 140 mm) Volume: 4 l	1 Piece (1 Box × 1 Piece)	E2396-6001
	 Create PDF Data Sheet		
	1 Litre Square Cooler Container, Green (270 x 195 x 115 mm) Volume: 1 l	1 Piece (1 Box × 1 Piece)	E2396-6002
	 Create PDF Data Sheet		
	9 Litre Oblong Cooler Container, Red (521 x 295 x 165 mm) Volume: 9 l	1 Piece (1 Box × 1 Piece)	E2396-6004
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk