



## Ultra High Recovery Microcentrifuge Tubes



### Product Highlights

- Max. centrifugation rate: 20,000 x g
- Temperature range: -80 °C to 121 °C
- Autoclavable
- Certified RNase, DNase, DNA and PCR Inhibitor free

## Product Applications

Sample preparation  
Centrifugation  
RNA or DNA pelleting  
Protein handling

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
Material	Polypropylene
DNA free	Yes
DNase free	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Natural
Operating temperature	-80°C to 121°C
g-Force max.	20,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Ultra High Recovery Microcentrifuge Tubes







Convenient for RNA or DNA pelleting or protein handling which often requires a low-adhesion surface, these ultra high recovery polypropylene tubes are completely non-reactive.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.5 ml Ultra High Recovery Microcentrifuge Tube</a> <b>Volume:</b> 0.5 ml  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1405-2600
 <a href="#">1.5 ml Ultra High Recovery Microcentrifuge Tube</a> <b>Volume:</b> 1.5 ml  <a href="#">Create PDF Data Sheet</a>	250 Tubes (1 Bag × 250 Tubes)	E1415-2600
 <a href="#">2.0 ml Ultra High Recovery Microcentrifuge Tube</a> <b>Volume:</b> 2.0 ml  <a href="#">Create PDF Data Sheet</a>	250 Tubes (1 Bag × 250 Tubes)	I1420-2600

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Cap-Locks for 1.5 / 2.0 ml Microcentrifuge Tubes</a>	100 Pcs. (1 Bag × 100 Pcs.)	I1415-1508

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)