



## 1.5 ml High Speed, Easy Close Microcentrifuge Tube



### Product Highlights

- Temperature range: -90 °C to 121 °C
- Autoclavable
- Certified RNase, DNase & DNA free

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Autoclavable	Yes
Volume	1.5 ml
Pack Size	500 Tubes (500 Tubes)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
RNase free	Yes
Color	Various
Operating temperature	-90 °C to 121 °C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 1.5 ml High Speed, Easy Close Microcentrifuge Tube













Solidly constructed for high speed centrifugation of max. 35,000 x g (not routinely tested at this force)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1.5 ml Copolymer Microcentrifuge Tube, Mixed</a> <b>Color:</b> Mixed  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-2508
 <a href="#">1.5 ml Copolymer Microcentrifuge Tube, Yellow</a> <b>Color:</b> Yellow  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-2506
 <a href="#">1.5 ml Copolymer Microcentrifuge Tube, Violet</a> <b>Color:</b> Violet  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-2505
 <a href="#">1.5 ml Copolymer Microcentrifuge Tube, Red</a> <b>Color:</b> Red  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-2504
 <a href="#">1.5 ml Copolymer Microcentrifuge Tube, Orange</a> <b>Color:</b> Orange  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-2503
 <a href="#">1.5 ml Copolymer Microcentrifuge Tube, Natural</a> <b>Color:</b> Natural  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-2500

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Cap-Locks for 1.5 / 2.0 ml Microcentrifuge Tubes</a>	100 Pcs. (1 Bag × 100 Pcs.)	I1415-1508

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)