




## **Alpha-Numeric Coded Tubes in Rack**



#### **Product Highlights**

- Refillable! Tube refills are supplied in a frame ready to load into racks
- The tube frame prevents tubes falling out of rack or the codes getting out of sequence
- SBS footprint compatible with all robotics applications
- Alpha-numeric codes printed on bottom of each tube
- Ideal for archiving samples
- Suitable for stackable systems
- Lockable, clear lid
- Temperature use -80 °C to 121 °C

### **Product Applications**

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







### **General Data**

Art. No.	See variations
----------	----------------

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







# More informations about Alpha-Numeric Coded Tubes in Rack

Alpha-numeric coded tubes in rack. Available with 0.65 ml or 1.2 ml tubes. Tube refills are supplied in a frame ready to load into racks. The tube frame prevents tubes falling out of rack or the codes getting out of sequence. Ideal for archiving samples and suitable for stackable systems.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD** 5 Tanners Drive

MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk







### **All Variations**

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Tube Rack with Lid, Pre-Filled with 96 x 0.65 ml Alpha-Numeric Coded Tubes  Packaging Unit: Racks  Volume: 0.65 ml  Create PDF Data Sheet	960 Tubes (10 Racks × 96 Tubes)	E1780-3100
	Tube Rack with Lid, Pre-Filled with 96 x 1.2 ml Alpha-Numeric Coded Tubes Packaging Unit: Racks Volume: 1.2 ml  Create PDF Data Sheet	960 Tubes (10 Racks × 96 Tubes)	E1780-2100
	Tube Rack Refill Frame, Loaded with 96 x 1.2 ml Alpha-Numeric Coded Tubes  Packaging Unit: Refills  Volume: 1.2 ml  Create PDF Data Sheet	960 Tubes (10 Refills × 96 Tubes)	E1780-2700

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.



5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk

