



Coloured Cap Inserts



Product Highlights

- Use to colour code the caps of Starlab's Cryogenic Vials
- Compatible with all Starlab Cryogenic Vials
- Offered in an array of colors for easy identification of samples

Product Applications

Colour coding for the caps of all Starlab Cryogenic Vials (except E3090-6321)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Sterile	No
Pack Size	500 Pcs. (1 Bag × 500 Pcs.)
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Various
Operating temperature	-196°C (only be used in the gas phase of liquid nitrogen)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Coloured Cap Inserts

For use with all Starlab Cryogenic Vials. Use to colour code your work and identify your samples more easily.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Coloured Cap Inserts for Cryogenic Vials, Blue Color: Blue  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6001
 Coloured Cap Inserts for Cryogenic Vials, Green Color: Green  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6002
 Coloured Cap Inserts for Cryogenic Vials, Orange Color: Orange  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6003
 Coloured Cap Inserts for Cryogenic Vials, Red Color: Red  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6004
 Coloured Cap Inserts for Cryogenic Vials, Violet Color: Violet  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6005
 Coloured Cap Inserts for Cryogenic Vials, Yellow Color: Yellow  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6006
 Coloured Cap Inserts for Cryogenic Vials, Brown Color: Brown  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6007

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Coloured Cap Inserts for Cryogenic Vials, White Color: White  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6009
 Coloured Cap Inserts for Cryogenic Vials, Mixed Color: Mixed  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6099
 Coloured Cap Inserts for Cryogenic Vials, Pink Color: Pink  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6104
 Coloured Cap Inserts for Cryogenic Vials, Grey Color: Gray  Create PDF Data Sheet	500 Pcs. (1 Bag × 500 Pcs.)	E1494-6108

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk