



## 0.2 ml PCR Tubes



### Product Highlights

- PCR Tubes manufactured from medical grade polypropylene
- Positive click on closing gives assurance that tubes are fully closed
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Certified free of RNase, Dnase, DNA & Pyrogen-free
- PCR inhibitor free
- Supplied in resealable bags

## Product Applications

Standard PCR (Endpoint PCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
Volume	0.2 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Various
Suitable for PCR	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 0.2 ml PCR Tubes











Starlab 0.2 ml [PCR tubes](#) with attached caps provide optimal sealing which minimises fluid loss through evaporation. The caps are easy to open and close yet ensure a tight fit to reduce sample evaporation.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)











## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.2 ml PCR Tube, Domed Cap, Natural</a> <b>Cap style:</b> Domed <b>Color:</b> Natural <b>Pack Size:</b> 1,000 Tubes (1 Bag × 1000 Tubes)  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (1 Bag × 1000 Tubes)	I1402-4300
 <a href="#">0.2 ml PCR Tube, Domed Cap, Mixed</a> <b>Cap style:</b> Domed <b>Color:</b> Mixed <b>Pack Size:</b> 1,000 Tubes (5 Bags × 200 Tubes)  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (5 Bags × 200 Tubes)	I1402-4308
 <a href="#">0.2 ml PCR Tube, Flat Cap, Natural</a> <b>Cap style:</b> Flat cap <b>Color:</b> Natural <b>Pack Size:</b> 1,000 Tubes (1 Bag × 1000 Tubes)  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (1 Bag × 1000 Tubes)	I1402-8100
 <a href="#">0.2 ml PCR Tube, Flat Cap, Mixed</a> <b>Cap style:</b> Flat cap <b>Color:</b> Mixed <b>Pack Size:</b> 1,000 Tubes (5 Bags × 200 Tubes)  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (5 Bags × 200 Tubes)	I1402-8108
 <a href="#">0.2 ml PCR Tube, Flat Cap (Xtra-Clear), Natural</a> <b>Cap style:</b> Flat (Xtra-Clear) <b>Color:</b> Natural <b>Pack Size:</b> 1,000 Tubes (1 Bag × 1000 Tubes)  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (1 Bag × 1000 Tubes)	I1402-8200

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.2 ml PCR Tube, Flat Cap, Natural (Case Size)</a> <b>Cap style:</b> Flat cap <b>Color:</b> Natural <b>Pack Size:</b> 10,000 Tubes (10 Bags × 1000 Tubes)  <a href="#">Create PDF Data Sheet</a>	10,000 Tubes (10 Bags × 1000 Tubes)	I1402-8100-C
 <a href="#">0.2 ml PCR Tube, Flat Cap, Mixed (Case Size)</a> <b>Cap style:</b> Flat cap <b>Color:</b> Mixed <b>Pack Size:</b> 10,000 Tubes (50 Bags × 200 Tubes)  <a href="#">Create PDF Data Sheet</a>	10,000 Tubes (50 Bags × 200 Tubes)	I1402-8108-C
 <a href="#">0.2ml PCR Tube, Domed Shielded Cap, Natural</a> <b>Cap style:</b> Domed shielded cap <b>Color:</b> Natural <b>Pack Size:</b> 1,000 Tubes (1 Bag × 1000 Tubes)  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (1 Bag × 1000 Tubes)	S1602-4300
 <a href="#">0.2 ml PCR Tube, Flat Shielded Cap, Natural</a> <b>Cap style:</b> Flat shielded cap <b>Color:</b> Natural <b>Pack Size:</b> 1,000 Tubes (1 Bag × 1000 Tubes)  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (1 Bag × 1000 Tubes)	B1402-5500

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)