



## PCR Strips with attached Caps



### Product Highlights

- Manufactured from medical grade polypropylene
- Attached caps eliminates the possibility of cross-contamination through accidental tube–cap mix up
- PCR inhibitor free
- Certified RNase, Dnase, DNA and Pyrogen-free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags

## Product Applications

Standard PCR (Endpoint PCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
No. of tubes	8
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Various
Suitable for PCR	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about PCR Strips with attached Caps

### Help eliminate cross-contamination with PCR Strips with Attached Caps

Discover convenience and efficiency with Starlab's PCR strips with attached caps, designed for mid-throughput settings. These PCR strips ensure all tubes are equally protected, supporting consistent results across samples, as well as simplifying handling to increase workflow efficiency.

With a selection of 0.1ml or 0.2ml tubes, all PCR strips feature ultra-thin, uniform walls to ensure optimal heat transfer. Choose between flat or domed caps, individually-attached or attached cap strip. For qPCR, select strips with Xtra-Clear caps. All PCR strips are certified free of RNase, DNase, DNA and PCR inhibitors, and tested Pyrogen free.




## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">0.2 ml 8-Strip PCR Tubes, Individually Attached Flat Caps (Xtra-Clear)</a>  <b>Cap style:</b> Individually attached flat caps (X-Clear)  <b>Volume:</b> 0.2 ml  <b>Pack Size:</b> 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)	A1402-3700
 <p><a href="#">0.2 ml 8-Strip PCR Tubes, Individually Attached Flat Caps (Xtra-Clear)(Case Size)</a>  <b>Cap style:</b> Individually attached flat caps (X-Clear)  <b>Volume:</b> 0.2 ml  <b>Pack Size:</b> 1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)	A1402-3700-C
 <p><a href="#">0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural</a>  <b>Cap style:</b> Individually attached domed caps  <b>Volume:</b> 0.2 ml  <b>Pack Size:</b> 120 Tube Strips with Caps (1 Packs × 120 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	120 Tube Strips with Caps (1 Packs × 120 Tube Strips with Caps)	I1402-2900
 <p><a href="#">8-Strip PCR Tubes with Individual Attached Caps, 0.2 ml (mixed)</a>  <b>Cap style:</b> Individually attached domed caps  <b>Volume:</b> 0.2 ml  <b>Pack Size:</b> 120 Tube Strips with Caps (5 Packs × 24 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	120 Tube Strips with Caps (5 Packs × 24 Tube Strips with Caps)	I1402-2908

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: infoline@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">0.1 ml 8-Strip PCR Tubes, Low-Profile, Individually Attached Flat Caps (Xtra-Clear)</a>  <b>Cap style:</b> Individually attached flat caps – (Xtra-Clear), low-profile  <b>Volume:</b> 0.1 ml  <b>Pack Size:</b> 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	<p>120 Tube Strips with Caps            (1 Pack × 120 Tube Strips with Caps)</p>	<p>A1402-3800</p>
 <p><a href="#">0.2 ml 8-Strip PCR Tubes, Individually Attached Domed Caps, Natural (Case Size)</a>  <b>Cap style:</b> Individually attached domed caps  <b>Volume:</b> 0.2 ml  <b>Pack Size:</b> 1,200 Tube Strips with Caps (10 Packs × 120 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	<p>1,200 Tube Strips with Caps            (10 Packs × 120 Tube Strips with Caps)</p>	<p>I1402-2900-C</p>
 <p><a href="#">0.2ml 8-Strip PCR Tubes, Individually Attached Domed Shielded Caps</a>  <b>Cap style:</b> Individually attached domed shielded caps  <b>Volume:</b> 0.2 ml  <b>Pack Size:</b> 120 Tube Strips with Caps (1 Pack × 120 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	<p>120 Tube Strips with Caps            (1 Pack × 120 Tube Strips with Caps)</p>	<p>S1602-2900</p>
 <p><a href="#">0.2 ml 8-Strip PCR Tubes, Natural with Attached Flat (Xtra-Clear) Cap Strip</a>  <b>Cap style:</b> Attached flat cap strip (Xtra-Clear)  <b>Volume:</b> 0.2 ml  <b>Pack Size:</b> 125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	<p>125 Tube Strips with Caps            (1 Pack × 125 Tube Strips with Caps)</p>	<p>I1402-2700</p>
 <p><a href="#">0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip</a>  <b>Cap style:</b> Attached domed cap strip  <b>Volume:</b> 0.2 ml  <b>Pack Size:</b> 125 Tube Strips with Caps (1 Pack × 125 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	<p>125 Tube Strips with Caps            (1 Pack × 125 Tube Strips with Caps)</p>	<p>I1402-2800</p>

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
 5 Tanners Drive  
 MK14 5BU Milton Keynes  
 Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">0.2 ml 8-Strip PCR Tubes, Natural with Attached Flat (Xtra-Clear) Cap Strip (Case Size)</a> <b>Cap style:</b> Attached flat cap strip (Xtra-Clear) <b>Volume:</b> 0.2 ml <b>Pack Size:</b> 1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)	I1402-2700-C
 <p><a href="#">0.2 ml 8-Strip PCR Tubes, Natural with Attached Domed Cap Strip (Case Size)</a> <b>Cap style:</b> Attached domed cap strip <b>Volume:</b> 0.2 ml <b>Pack Size:</b> 1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	1,250 Tube Strips with Caps (10 Packs × 125 Tube Strips with Caps)	I1402-2800-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)