



## Real-Time PCR Strips with attached Caps



### Product Highlights

- Manufactured from medical grade polypropylene
- 'Non-flex' PCR Strips for easier handling
- Xtra-Clear caps for qPCR
- Certified RNase, Dnase, DNA & Pyrogen-free
- PCR Inhibitor free
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- Supplied in resealable bags

## Product Applications

Standard PCR (Endpoint PCR)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
Pack Size	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Various
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Real-Time PCR Strips with attached Caps











Starlab non-flex tubes have a robust, triple link between each tube making them less flexible and easier to handle. The individually attached caps eliminates possible cross contamination. The caps are set an angle which allows the tubes to be nested together in racks without the lids interfering with one another.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">0.2 ml 8-Strip Non-Flex PCR Tubes, Natural, Individually Attached Flat Caps (Xtra-Clear)</a>  <b>Volume:</b> 0.2 ml  <b>Color:</b> Natural</p> <p> <a href="#">Create PDF Data Sheet</a></p>	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1402-3700
 <p><a href="#">0.1 ml 8-Strip Non-Flex PCR Tubes, Low-Profile, Natural, Individually Attached Flat Caps (Xtra-Clear)</a>  <b>Volume:</b> 0.1 ml  <b>Color:</b> Natural</p> <p> <a href="#">Create PDF Data Sheet</a></p>	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1402-3800
 <p><a href="#">0.2 ml 8-Strip Non-Flex PCR Tubes, White, Individually Attached Flat Caps (Xtra-Clear)</a>  <b>Volume:</b> 0.2 ml  <b>Color:</b> White</p> <p> <a href="#">Create PDF Data Sheet</a></p>	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1409-3700
 <p><a href="#">0.1 ml 8-Strip Non-Flex PCR Tubes, Low-Profile, White, Individually Attached Flat Caps (Xtra-Clear)</a>  <b>Volume:</b> 0.1 ml  <b>Color:</b> White</p> <p> <a href="#">Create PDF Data Sheet</a></p>	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1409-3800
 <p><a href="#">0.2 ml 8-Strip Non-Flex PCR Tubes, Mixed, Individually Attached Flat Caps (Xtra-Clear)</a>  <b>Volume:</b> 0.2 ml  <b>Color:</b> Mixed</p> <p> <a href="#">Create PDF Data Sheet</a></p>	120 Tube Strips with Caps (1 Bag × 120 Tube Strips with Caps)	I1402-3708

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)