



PCR Strip Tube & Cap Combo Packs



Product Highlights

- Manufactured from medical grade polypropylene
- Strip PCR tube and PCR strip cap combo pack
- Endotoxin free.
- Supplied in resealable bags

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Sterile	No
Autoclavable	No
Pack Size	125 Pairs (1 × 125 Pairs)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
RNase free	Yes
Suitable for PCR	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk




More informations about PCR Strip Tube & Cap Combo Packs

Buy your PCR strip tubes and strip caps together! Supplied unassembled, the tube strips and cap strips are packaged separately in resealable bags.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.1 ml 8-Strip PCR Tubes, Low-Profile, Natural and 8-Strip PCR Caps, Flat (Xtra-Clear) Volume: 0.1 ml Cap style: Flat (Xtra-Clear) Color: Natural  Create PDF Data Sheet	125 Pairs (1 Pack × 125 Pairs)	I1402-4700
 0.2 ml 8-Strip PCR Tubes, Natural and 8-Strip PCR Caps, Domed Volume: 0.2 ml Cap style: Domed Color: Natural  Create PDF Data Sheet	125 Pairs (1 Pack × 125 Pairs)	I1402-6700
 0.2 ml 8-Strip PCR Tubes, Natural and 8-Strip PCR Caps, Flat (Xtra-Clear) Volume: 0.2 ml Cap style: Flat (Xtra-Clear) Color: Natural  Create PDF Data Sheet	125 Pairs (1 Pack × 125 Pairs)	I1402-5700
 0.1 ml 8-Strip PCR Tubes, Low-Profile, White and 8-Strip PCR Caps, Flat (Xtra-Clear) Volume: 0.1 ml Cap style: Flat (Xtra-Clear) Color: White  Create PDF Data Sheet	125 Pairs (1 Box × 125 Pairs)	I1409-4700
 0.2 ml 8-Strip PCR Tubes, White and 8-Strip PCR Caps, Flat (Xtra-Clear) Volume: 0.2 ml Cap style: Flat (Xtra-Clear) Color: White  Create PDF Data Sheet	125 Pairs (1 Box × 125 Pairs)	I1409-5700

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk