



Deepwell Plates



Product Highlights

- Virgin polypropylene plates with SBS footprint
- Alpha-numeric reference grid (moulded)
- Plate temperature use: -80 °C to 121 °C.
- Autoclavable
- Max. centrifugation rate: 5000 x g

Product Applications

Sample preparation and storage
DNA and RNA isolation with high throughput
Plasmid purification with high throughput
Storage of cDNA or genomic banks
Storage of siRNA or oligonucleotide libraries
Protein analysis
Bacteria and yeast cultivation
Dilutions series

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Autoclavable	Yes
Material	Polypropylene
DNA free	Yes
DNase free	Yes
RNase free	Yes
g-Force max.	5,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Deepwell Plates











Virgin polypropylene plates with SBS footprint. Alpha-numeric reference grid (moulded).

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk















All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 2.5ml 96 Deepwell Plate, Square Wells with Round Bottoms Volume: 2.5 ml Sterile: No Pack Size: 20 Plates (1 Box × 20 Plates)  Create PDF Data Sheet	20 Plates (1 Box × 20 Plates)	E2896-2500
 10 ml 24 Deepwell Plate, Square Wells with V-Shape Bottoms Volume: 10 ml Sterile: No Pack Size: 20 Plates (1 Box × 20 Plates)  Create PDF Data Sheet	20 Plates (1 Box × 20 Plates)	E2824-0000
 5.0 ml 48 Deepwell Plate, Rectangle Wells with V-Shape Bottoms Volume: 5.0 ml Sterile: No Pack Size: 20 Plates (1 Box × 20 Plates)  Create PDF Data Sheet	20 Plates (1 Box × 20 Plates)	E2848-0000
 1.2 ml 96 Deepwell Plate, Round Wells with Round Bottoms, Low-Profile Volume: 1.2 ml Sterile: No Pack Size: 50 Plates (1 Box × 50 Plates)  Create PDF Data Sheet	50 Plates (1 Box × 50 Plates)	E2896-0120
 2.2 ml 96 Deepwell Plate, Square Wells with V-Shaped Bottoms Deepwell Plates Volume: 2.20 ml Sterile: No Pack Size: 50 Plates (1 Box × 50 Plates)  Create PDF Data Sheet	50 Plates (1 Box × 50 Plates)	E2896-0220

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 2.0 ml 96 Deepwell Plate, Round Wells with Round Bottoms Volume: 2.0 ml Sterile: No Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E2896-2100
 1.2 ml 96 Deepwell Plate, Square Wells with Conical Bottoms Volume: 1.2 ml Sterile: No Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E2896-1820
 2.2 ml 96 Deepwell Plate, Square Wells with Conical Bottoms Volume: 2.2 ml Sterile: No Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E2896-1810
 0.6 ml 96 Deepwell Plate, Round Wells with Conical Bottoms Volume: 0.6 ml Sterile: No Pack Size: 5 Plates (1 Box × 5 Plates)  Create PDF Data Sheet	5 Plates (1 Box × 5 Plates)	E2896-0600
 1.1 ml 96 Deepwell Plate, Round Wells with Round Bottoms, Low-Profile (Sterile) Volume: 1.1 ml Sterile: Yes Pack Size: 10 Plates (2 Packs × 5 Plates)  Create PDF Data Sheet	10 Plates (2 Packs × 5 Plates)	E2896-0110
 240 µl 384 Deepwell Plate, Square Wells with V-Shape Bottoms Volume: 0.24 ml Sterile: No Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet	10 Plates (1 Box × 10 Plates)	E2384-2400

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk