



200 µl XL Filter Tip, sterile, Ind-Wrap/ 100 tips (100 x individually wrapped sachets)



Product Highlights

- Compatible with a wide range of pipettes
- Made from polypropylene
- For some tips sterile racks or bags available

Product Applications

Pipetting liquids from and into long tubes
Phase extraction and removal of lower phases

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	I1011-0070
Pack Size	100 Individually wrapped sachets (1 Bag × 100 Individually wrapped sachets)
Sterile	Yes
Filter tip	Filter
Volume range	5 - 200 µl
Volume	200 µl
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk




More informations about 200 µl XL Filter Tip, sterile, Ind-Wrap/ 100 tips (100 x individually wrapped sachets)

Extra long tips for easy handling of liquids in long tubes. Phase extraction and removal of lower phases. 1 pack of 100 individually-wrapped tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1.25 ml XL Graduated Filter Tip, Racked (sterile) Volume: 1.25 ml Sterile: Yes Pack Size: 768 Tips (8 Racks × 96 Tips)  Create PDF Data Sheet	768 Tips (8 Racks × 96 Tips)	I1011-3010
 200 µl XL Filter Tip, sterile, Racked Volume: 200 µl Sterile: Yes Pack Size: 800 Tips (8 Racks × 100 Tips)  Create PDF Data Sheet	800 Tips (8 Racks × 100 Tips)	I1011-5830
 200 µl XL Filter Tip, sterile, Racked Volume: 200 µl Sterile: Yes Pack Size: 1,632 Tips (8 Racks × 204 Tips)  Create PDF Data Sheet	1,632 Tips (8 Racks × 204 Tips)	I1011-0830
 1,000 µl XL Bevelled Filter Tip, Loose / 1000 tips (1 bag x 1000 tips) Volume: 1000 µl Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	E1011-1008
 200 µl XL Filter Tip, sterile, Ind-Wrap/ 100 tips (100 x individually wrapped sachets) Volume: 200 µl Sterile: Yes Pack Size: 100 Individually wrapped sachets (1 Bag × 100 Individually wrapped sachets)  Create PDF Data Sheet	100 Individually wrapped sachets (1 Bag × 100 Individually wrapped sachets)	I1011-0070

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk