



## Disposable PS Reservoirs (StarTub PS)



### Product Highlights


- White polystyrene (PS) reservoir
- Sterile and Non-Sterile
- For Single use
- Capacity: 25 - 100 ml

## Product Applications

Reagent takeup with multichannel pipette

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	55 mm × 65 mm × 26 mm
Dimensions (W × D)	55 mm × 65 mm
Autoclavable	No
Material	Polystyrene
Color	White
Bottoms	Conical
Reusable	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Disposable PS Reservoirs (StarTub PS)

55 ml polystyrene, sterile individually wrapped reservoirs with V-shaped bottom.









Dimensions (WxDxH): 55ml: 149 x 65 x 26 mm

100 ml polystyrene bulk-packed reservoirs with trough-shaped bottom.

Dimensions (WxDxH): 144 x 85 x 33 mm



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">StarTub PS Reagent Reservoir, 55 ml, Individually Wrapped (Sterile)</a> <b>Sterile:</b> Yes <b>Volume:</b> 55 ml   <a href="#">Create PDF Data Sheet</a>	50 Pcs. (1 Box × 50 Pcs.)	E2310-1010
	<a href="#">StarTub PS Reagent Reservoir, 100 ml, Bulk</a> <b>Sterile:</b> No <b>Volume:</b> 100 ml   <a href="#">Create PDF Data Sheet</a>	100 Pcs. (1 Box × 100 Pcs.)	E2310-2000
	<a href="#">StarTub PS Reagent Reservoir, 25 ml, Individually Wrapped (Sterile)</a> <b>Sterile:</b> Yes <b>Volume:</b> 25 ml   <a href="#">Create PDF Data Sheet</a>	50 Pcs. (1 Box × 50 Pcs.)	E2310-1025
	<a href="#">StarTub PS Reagent Reservoir, 100 ml, Individually Wrapped (Sterile)</a> <b>Sterile:</b> Yes <b>Volume:</b> 100 ml   <a href="#">Create PDF Data Sheet</a>	50 Pcs. (1 Box × 50 Pcs.)	E2310-1100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)