



ChannelMate Reservoirs



Product Highlights

- Non-tipping, durable base features high-contrast volume graduations
- Skirt-less polystyrene inserts require up to 40% less plastic and create minimal waste
- Nested inserts save up to 65% of the storage space required with other reservoirs
- Two systems: 25 ml & 100 ml (bases and inserts can be purchased separately)

Product Applications

Reagent takeup with multichannel pipette

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	See variations
Sterile	Yes
Material	Polystyrene
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk




More informations about ChannelMate Reservoirs













Reagent reservoir refills save waste and storage space. ChannelMate™ combines a sturdy reusable base with ultra-clear inserts to create a unique refillable reservoir system. DNase, RNase, and pyrogen free. Sterile.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 ChannelMate Insert, 25 ml (Sterile) Component Type: Insert Volume: 25 ml  Create PDF Data Sheet	50 Pcs. (1 Box × 50 Pcs.)	E1346-2510
 ChannelMate Base, 25 ml (Sterile) Component Type: Base Volume: 25 ml  Create PDF Data Sheet	10 Pcs. (1 Box × 10 Pcs.)	E1306-2590
 ChannelMate System, 25 ml (Sterile) Component Type: System Volume: 25 ml  Create PDF Data Sheet	200 Pcs. (1 Box × 200 Pcs.)	E1306-2510
 ChannelMate Insert, 100 ml (Sterile) Component Type: Insert Volume: 100 ml  Create PDF Data Sheet	50 Pcs. (1 Box × 50 Pcs.)	E1346-1010
 ChannelMate Base, 100 ml (Sterile) Component Type: Base Volume: 100 ml  Create PDF Data Sheet	10 Pcs. (1 Box × 10 Pcs.)	E1306-1090
 ChannelMate System, 100 ml (Sterile) Component Type: System Volume: 100 ml  Create PDF Data Sheet	200 Pcs. (1 Box × 200 Pcs.)	E1306-1010

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk