




Standard cuvettes



Product Highlights

- Polystyrene
- Cuvettes supplied in polystyrene rack with lid for safe storage.

Product Applications

For assays in the visible spectral range 340 to 750 nm (E1412-4160 range 340 to 900 nm).

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Material	Polystyrene

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk





More informations about Standard cuvettes

Disposable [cuvettes](#) suitable for use with most open spectrophotometers. Provides optimum transmission level on the whole visible spectral. The material used avoids any possible measurement interference. High reliability is assured, with maximum absorption variation at: $\pm 1\%$.

Light beam width: for 10 mm. The two sides not crossed by the light beam are ribbed to an easy identification of the [cuvette](#) position inside the spectrophotometer measurement chamber, resulting in an easy positioning and removal.



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1.5 ml SemiMicro Cuvette, Square, Polystyrene Volume: 1.5 ml  Create PDF Data Sheet	500 Pcs. (5 Packs × 100 Pcs.)	E1412-4150

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk