



Screw Cap Containers, non-sterile



Product Highlights

Product Applications

Sample preparation

Sample storage

Preparation of master mixes and buffers

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|--------------|----------------|
| Art. No. | See variations |
| Sterile | No |
| Autoclavable | No |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Screw Cap Containers, non-sterile





60 ml crystal clear polystyrene container with red polyethylene screw cap. Aseptically* produced. Temperature use: -10 °C to 70 °C.

125 ml polypropylene container with red polyethylene screw cap. 20 ml graduations marked up to 100 ml. Frosted writing panel. Aseptically* produced. Temperature use: -10 °C to 70 °C.

** Aseptically produced in a fully-automated and bio-controlled environment.*



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|-----------------------------|------------|
|  60 ml Polystyrene Screw Cap Container, non-sterile Volume: 60 ml  Create PDF Data Sheet | 600 Pcs. (1 Box × 600 Pcs.) | E1412-6004 |
|  125 ml Polypropylene Graduated Screw Cap Container, non-sterile Volume: 125 ml  Create PDF Data Sheet | 475 Pcs. (1 Box × 475 Pcs.) | E1412-6104 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk