



## CytoOne® Dishes



### Product Highlights

- Sturdy, flat bases resist bowing and warping
- Optical clarity provides distortion-free microscopy
- Stacking rings with vented design for thermal equalisation
- Vented lids improve gas exchange
- Gamma sterilised (SAL  $10^{-5}$ ) and certified RNase, DNase, DNA and pyrogen free
- All CytoOne® dishes have untreated lids to minimise condensation
- Easy to open sleeves in convenient, small pack sizes

## Product Applications

Expansion and culture of adherent and suspension cells  
Performance of cell-based assays and microscopic analysis

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Sterile	Yes
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Clear
Non-cytotoxic	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about CytoOne® Dishes















### Secure handling

Available with [tissue culture](#) treatment for optimum cell attachment and growth or non-treated for suspension cultures. Bevelled outer edges on the 35 mm and 60 mm dish bases make these small dishes easier to handle without accidentally displacing the lids, even when stacking. The lids of the 150 mm dishes have four ergonomically placed grip tags for secure handling with less strain.

- › Stacking rings with vented design for thermal equalisation
- › Vented lids improve gas exchange
- › All [CytoOne](#)® dishes have untreated lids to minimise condensation
- › Easy to open sleeves in convenient, small pack sizes





## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">150 x 20 mm CytoOne® Dish, TC-Treated</a> <b>Working volume:</b> 30 ml <b>Treatment:</b> TC Treated   <a href="#">Create PDF Data Sheet</a>	60 Dishes (12 Packets × 5 Dishes)	CC7682-3614
 <a href="#">100 x 20 mm CytoOne® Dish, TC-Treated</a> <b>Working volume:</b> 10 ml <b>Treatment:</b> TC Treated   <a href="#">Create PDF Data Sheet</a>	300 Dishes (30 Packets × 10 Dishes)	CC7682-3394
 <a href="#">60 x 15 mm CytoOne® Dish, TC-Treated</a> <b>Working volume:</b> 5.0 ml <b>Treatment:</b> TC Treated   <a href="#">Create PDF Data Sheet</a>	300 Dishes (30 Packets × 10 Dishes)	CC7682-3359
 <a href="#">35 x 10 mm CytoOne® Dish, TC-Treated</a> <b>Working volume:</b> 2.0 ml <b>Treatment:</b> TC Treated   <a href="#">Create PDF Data Sheet</a>	300 Dishes (30 Packets × 10 Dishes)	CC7682-3340
 <a href="#">150 x 20 mm CytoOne® Dish, Non-Treated</a> <b>Working volume:</b> 30 ml <b>Treatment:</b> Non-Treated   <a href="#">Create PDF Data Sheet</a>	60 Dishes (12 Packets × 5 Dishes)	CC7672-3614
 <a href="#">100 x 20 mm CytoOne® Dish, Non-Treated</a> <b>Working volume:</b> 10 ml <b>Treatment:</b> Non-Treated   <a href="#">Create PDF Data Sheet</a>	300 Dishes (30 Packets × 10 Dishes)	CC7672-3394
 <a href="#">60 x 15 mm CytoOne® Dish, Non-Treated</a> <b>Working volume:</b> 5.0 ml <b>Treatment:</b> Non-Treated   <a href="#">Create PDF Data Sheet</a>	300 Dishes (30 Packets × 10 Dishes)	CC7672-3359

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">35 x 10 mm CytoOne® Dish, Non-Treated</a> <b>Working volume:</b> 2.0 ml <b>Treatment:</b> Non-Treated</p> <p> <a href="#">Create PDF Data Sheet</a></p>	300 Dishes (30 Packets × 10 Dishes)	CC7672-3340

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)