



CytoOne® Flasks with 2-Position Plug Seal Caps



Product Highlights

- All flasks 100% pressure tested for leak-free assurance
- Stacking rims with vented design for thermal equalisation
- Anti-tilt skirts and base bars improve stability
- Wide necks provide easy access while minimising contamination
- Gamma sterilised (SAL 10^{-5}) and certified RNase, DNase, DNA and Pyrogen free
- Optical clarity provides distortion-free microscopy
- Easy open, resealable zip top bags protect unused flasks

Product Applications

Expansion and culture of adherent and suspension cells
Performance of cell-based assays and microscopic analysis

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Sterile	Yes
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Non-cytotoxic	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about CytoOne® Flasks with 2-Position Plug Seal Caps

Pressure tested















Choose tissue culture treated for optimum cell attachment and growth or non-treated for hybridoma and suspension cultures. The two-position plug seal cap lets you choose secure manual venting or an airtight seal. Flasks feature convenient moulded volume graduations and frosted writing surfaces. The T-25, T-75, and T-150 flasks have traceable lot numbers etched into the base.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk





All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 T-25 CytoOne® Flask, TC-Treated Flask Type: T-25 Treatment: TC Treated  Create PDF Data Sheet	300 Flasks (30 Bags × 10 Flasks)	CC7682-4325
 T-75 CytoOne® Flask, TC-Treated Flask Type: T-75 Treatment: TC Treated  Create PDF Data Sheet	100 Flasks (20 Bags × 5 Flasks)	CC7682-4175
 T-150 CytoOne® Flask, TC-Treated Flask Type: T-150 Treatment: TC Treated  Create PDF Data Sheet	40 Flasks (8 Bags × 5 Flasks)	CC7682-4415
 T-225 CytoOne® Flask, TC-Treated Flask Type: T-225 Treatment: TC Treated  Create PDF Data Sheet	25 Flasks (5 Bags × 5 Flasks)	CC7682-4225
 T-25 CytoOne® Flask, Non-Treated Flask Type: T-25 Treatment: Non-Treated  Create PDF Data Sheet	300 Flasks (30 Bags × 10 Flasks)	CC7672-4325
 T-75 CytoOne® Flask, Non-Treated Flask Type: T-75 Treatment: Non-Treated  Create PDF Data Sheet	100 Flasks (20 Bags × 5 Flasks)	CC7672-4175
 T-150 CytoOne® Flask, Non-Treated Flask Type: T-150 Treatment: Non-Treated  Create PDF Data Sheet	40 Flasks (8 Bags × 5 Flasks)	CC7672-4415

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>T-225 CytoOne® Flask, Non-Treated Flask Type: T-225 Treatment: Non-Treated</p> <p> Create PDF Data Sheet</p>	25 Flasks (5 Bags × 5 Flasks)	CC7672-4225

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk