



15 ml Centrifuge Tubes



Product Highlights

- Clear, sterile polypropylene tubes with green HDPE screw caps
- White writing surface and clear, printed graduations
- Perfect for centrifugation of cells: USP Class VI, non-cytotoxic and sterile (SAL 10^{-6}) material
- Certified RNase, DNase, DNA, heavy-metal and Pyrogen free

Product Applications

Centrifugation
Sample preparation
Preparation of master mixes and buffers
Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|----------------------|------------------|
| Art. No. | See variations |
| Dimensions (W × H) | 22 mm × 120.5 mm |
| Sterile | Yes |
| Autoclavable | Yes |
| Volume | 15 ml |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| Free of endotoxins | Yes |
| Free of heavy metals | Yes |
| Human DNA free | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Color | Clear |
| g-Force max. | 12,000 x g |
| Non-cytotoxic | Yes |
| Tube Base | Conical |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk







More informations about 15 ml Centrifuge Tubes

STARLAB's new [centrifuge tubes](#) are designed for the centrifugation of cells, but can also be used for standard centrifugation, preparation and storage applications. USP Class VI polymers and certified non-toxicity ensure the viability and integrity of cells. The sterility of assurance level (SAL) of 10^{-6} fulfills highest demands on sterility, meeting EN 556-1 standard (Sterilisation of medical devices). Racked tubes are delivered in recyclable cardboard racks.

Maximum centrifugation rates: Conical Tubes 15 & 50 ml: 12,000 x g and Skirted Tube 50 ml: 6,000 x g



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------------|------------|
|  15 ml Centrifuge Tube, Conical, Loose (Sterile) Packaging Unit: Bags  Create PDF Data Sheet | 500 Tubes (20 Bags × 25 Tubes) | E1415-0100 |
|  15 ml Centrifuge Tube, Conical, Racked (Sterile) Packaging Unit: Racks  Create PDF Data Sheet | 500 Tubes (20 Racks × 25 Tubes) | E1415-0700 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk