



## Autoclave bags



### Product Highlights

- Polypropylene bags for autoclaving up to 136°C
- 48 microns
- Printed with the Biohazard Symbol
- Use for contaminated waste disposal
- White write-on panel

## Product Applications

Disposal

Occupational health and safety

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Material	Polypropylene
Pack Size	200 Pcs. (200 Pcs.)
Autoclavable	Yes
Color	Clear

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Autoclave bags







Autoclave bags

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">Autoclave bags 310 x 660 mm</a> <b>Dimensions (W x H):</b> 310 mm x 660 mm	200 Pcs. (1 Box x 200 Pcs.)	E3300-3166
	<a href="#">Create PDF Data Sheet</a>		
	<a href="#">Autoclave bags 415 x 600 mm</a> <b>Dimensions (W x H):</b> 415 mm x 600 mm	200 Pcs. (1 Box x 200 Pcs.)	E3300-4070
	<a href="#">Create PDF Data Sheet</a>		
	<a href="#">Autoclave bags 610 x 810 mm</a> <b>Dimensions (W x H):</b> 610 mm x 810 mm	200 Pcs. (1 Box x 200 Pcs.)	E3300-6078
	<a href="#">Create PDF Data Sheet</a>		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)