



Autoclave bags



Product Highlights

- Polypropylene bags for autoclaving up to 136°C
- 48 microns
- Printed with the Biohazard Symbol
- Use for contaminated waste disposal
- White write-on panel

Product Applications

Disposal

Occupational health and safety

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Autoclavable	Yes
Pack Size	200 Pcs. (200 Pcs.)
Material	Polypropylene
Color	Clear

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Autoclave bags







Autoclave bags

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Autoclave bags 310 x 660 mm Dimensions (W x H): 310 mm x 660 mm  Create PDF Data Sheet	200 Pcs. (1 Box x 200 Pcs.)	E3300-3166
 Autoclave bags 415 x 600 mm Dimensions (W x H): 415 mm x 600 mm  Create PDF Data Sheet	200 Pcs. (1 Box x 200 Pcs.)	E3300-4070
 Autoclave bags 610 x 810 mm Dimensions (W x H): 610 mm x 810 mm  Create PDF Data Sheet	200 Pcs. (1 Box x 200 Pcs.)	E3300-6078

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk