




Round Gel-Loading Filter Tips



Product Highlights

- Flexible, long and narrow
- Loading samples onto gels (eg. polyacrylamide or agarose)
- Handling of small volumes
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Filter tip	Filter
Tip Type	round gel-loading

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Round Gel-Loading Filter Tips









All tips are certified RNase and DNase free. Packaging of these tips differs from the [TipOne](#)[®] range.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 10 µl Round Gel-Loading Filter Tip Volume: 10 µl Sterile: Yes  Create PDF Data Sheet	1,020 Tips (5 Racks × 204 Tips)	I1010-3800
 20 µl Round Gel-Loading Filter Tip (Sterile) Volume: 20 µl Sterile: Yes  Create PDF Data Sheet	1,020 Tips (5 Racks × 204 Tips)	E1021-0618
 100 µl Round Gel-Loading Filter Tip Volume: 100 µl Sterile: No  Create PDF Data Sheet	1,000 Tips (1 Bag × 1000 Tips)	E1022-0008
 100 µl Round Gel-Loading Filter Tip (Sterile) Volume: 100 µl Sterile: Yes  Create PDF Data Sheet	1,020 Tips (5 Racks × 204 Tips)	I1022-0810

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk