




Extended Length Tips



Product Highlights

- Compatible with a wide range of pipettes
- Made from polypropylene
- For some tips sterile racks or bags available

Product Applications

Pipetting liquids from and into long tubes
Phase extraction and removal of lower phases

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|------------|----------------|
| Art. No. | See variations |
| Filter tip | Non-filter |
| Material | Polypropylene |
| Color | Natural |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Extended Length Tips











Extra long tips for easy handling of liquids in long tubes. Phase extraction and removal of lower phases.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|----------------------------------|------------|
|  200 µl XL Tip Racked (sterile) Volume: 200 µl Sterile: Yes Packaging Unit: Racks  Create PDF Data Sheet | 1,632 Tips (8 Racks × 204 Tips) | E1011-0921 |
|  1,000 µl XL Bevelled Tip, Loose (non-sterile) Volume: 1000 µl Sterile: No Packaging Unit: Bag  Create PDF Data Sheet | 1,000 Tips (1 Bag × 1000 Tips) | E1011-1002 |
|  200 µl XL Tip, Loose (non-sterile) Volume: 200 µl Sterile: No Packaging Unit: Bag  Create PDF Data Sheet | 1,000 Tips (1 Bag × 1000 Tips) | I1011-0020 |
|  200 µl XL Tip, Racked (non-sterile) Volume: 200 µl Sterile: No Packaging Unit: Racks  Create PDF Data Sheet | 1,632 Tips (8 Racks × 204 Tips) | I1011-0920 |
|  1,000 µl XL Bevelled Tip, Racked (non-sterile) Volume: 1000 µl Sterile: No Packaging Unit: Racks  Create PDF Data Sheet | 1,000 Tips (10 Racks × 100 Tips) | E1011-1430 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk