



## 200 µl Bevelled TipOne® Tip, Natural



### Product Highlights

- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Certified RNase, DNase, DNA and Pyrogen free.

## Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 200 µl
Volume	200 µl
Tip Type	bevelled
Material	Polypropylene
Graduations	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 200 µl Bevelled TipOne® Tip, Natural

### Experience effortless precision with [TipOne® Pipette Tips](#)

Elevate your daily lab work with [TipOne® pipette tips](#)! Enjoy precise, reliable pipetting every time—no more wasted samples or frustrating fit issues. [TipOne®](#) tips are compatible with all popular pipettes for accurate results you can trust, every single day.

### Why [TipOne®](#) makes labs better:

**Universal tip fit:** Simplify your bench and inventory. The tip collar of [TipOne®](#) tips is engineered to fit all popular pipettes without any loss of accuracy or precision.

**Made for precision:** Expert design and the use of high-quality virgin materials means our tips provide consistent accuracy to help you get reliable results and reduce costly rework.

**Pure quality:** Sterile [TipOne®](#) tips are validated according to EN ISO 11137-2:2007. All [TipOne®](#) tips are certified RNase, DNase, DNA and pyrogen free.

### Flexible Packaging:

**Racked:** Sturdy, reusable racks stay stable on your bench. Quickly reload with [TipOne®](#) refills for reduced plastic waste and cost. Sterile options available.

**Refill:** Load tips into your racks with a few easy steps. Refills are available sterile or non-sterile, helping you cut down on plastic waste (and expenses).

**Stack Rack:** Handle large experiments effortlessly—dispense 960 tips from one compact unit that is perfect for multi-channel pipettes and busy workflows. Sterile options available.

**Loose:** For teaching or routine work, stock up on affordable quality tips. Supplied in convenient resealable bags.













### Join our journey to sustainable pipetting!

Reduce your lab's plastic waste! [TipOne®](#) racks and refill system components use recycled materials and are fully recyclable (except shrink wrap, loose bag, and tips) via the [TipOne®](#) recycling service.

We'll collect your used [TipOne®](#) plastic, recycle it and use the recycled plastic to make more [TipOne®](#) products. No more guessing if or how your recycling gets reused—[TipOne®](#) closes the loop for you!















## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (non-sterile)</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Racks × 96 Tips)	S1111-1800
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (non-sterile), Case</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	7,680 Tips (80 Racks × 96 Tips)	S1111-1800-C
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (sterile)</a>  <b>Sterile:</b> Yes  <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Racks × 96 Tips)	S1111-1810
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (sterile), Case</a>  <b>Sterile:</b> Yes  <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	7,680 Tips (80 Racks × 96 Tips)	S1111-1810-C
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (non-sterile)</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (10 Refills × 96 Tips)	S1111-1700
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (non-sterile), Case</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 9,600 Tips (100 Refills × 96 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	9,600 Tips (100 Refills × 96 Tips)	S1111-1700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)







PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack rack (non-sterile)</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 960 Tips (1 Stack rack × 960 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (1 Stack rack × 960 Tips)	S1111-1200
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack racks (non-sterile), Case</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 9,600 Tips (10 Stack racks × 960 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	9,600 Tips (10 Stack racks × 960 Tips)	S1111-1200-C
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack rack (sterile)</a>  <b>Sterile:</b> Yes  <b>Pack Size:</b> 960 Tips (1 Stack rack × 960 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	960 Tips (1 Stack rack × 960 Tips)	S1111-1210
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack racks (sterile), Case</a>  <b>Sterile:</b> Yes  <b>Pack Size:</b> 9,600 Tips (10 Stack racks × 960 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	9,600 Tips (10 Stack racks × 960 Tips)	S1111-1210-C
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Bag (non-sterile)</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 1,000 Tips (1 Bag × 1000 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	1,000 Tips (1 Bag × 1000 Tips)	S1111-1000
 <p><a href="#">200 µl Bevelled TipOne® Tip, Natural, Bags (non-sterile), Case</a>  <b>Sterile:</b> No  <b>Pack Size:</b> 25,000 Tips (25 Bags × 1000 Tips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	25,000 Tips (25 Bags × 1000 Tips)	S1111-1000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1111-1710
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 7,680 Tips (80 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Refills × 96 Tips)	S1111-1710-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)