



200 µl UltraPoint® graduated TipOne® Tip, Yellow



Product Highlights

- **Graduated** at 10 µl, 50 µl and 100 µl
- UltraPoint® thin-wall orifice helps prevent sample hang up
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Certified RNase, DNase, DNA and Pyrogen free.

Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|--------------|----------------|
| Art. No. | See variations |
| Autoclavable | Yes |
| Filter tip | Non-filter |
| Volume range | 1.0 - 200 µl |
| Volume | 200 µl |
| Material | Polypropylene |
| Graduations | Yes |
| Color | Yellow |
| Tip Type | UltraPoint® |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 200 µl UltraPoint® graduated TipOne® Tip, Yellow

Experience effortless precision with TipOne® Pipette Tips

Elevate your daily lab work with TipOne® pipette tips! Enjoy precise, reliable pipetting every time—no more wasted samples or frustrating fit issues. TipOne® tips are compatible with all popular pipettes for accurate results you can trust, every single day.

Why TipOne® makes labs better:

Universal tip fit: Simplify your bench and inventory. The tip collar of TipOne® tips is engineered to fit all popular pipettes without any loss of accuracy or precision.

Made for precision: Expert design and the use of high-quality virgin materials means our tips provide consistent accuracy to help you get reliable results and reduce costly rework.

Pure quality: Sterile TipOne® tips are validated according to EN ISO 11137-2:2007. All TipOne® tips are certified RNase, DNase, DNA and pyrogen free.

Flexible Packaging:

Racked: Sturdy, reusable racks stay stable on your bench. Quickly reload with TipOne® refills for reduced plastic waste and cost. Sterile options available.

Refill: Load tips into your racks with a few easy steps. Refills are available sterile or non-sterile, helping you cut down on plastic waste (and expenses).

Stack Rack: Handle large experiments effortlessly—dispense 960 tips from one compact unit that is perfect for multi-channel pipettes and busy workflows. Sterile options available.

Loose: For teaching or routine work, stock up on affordable quality tips. Supplied in convenient resealable bags.













Join our journey to sustainable pipetting!

Reduce your lab's plastic waste! TipOne® racks and refill system components use recycled materials and are fully recyclable (except shrink wrap, loose bag, and tips) via the TipOne® recycling service.

We'll collect your used TipOne® plastic, recycle it and use the recycled plastic to make more TipOne® products. No more guessing if or how your recycling gets reused—TipOne® closes the loop for you!















All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|------------------------------------|--------------|
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (non-sterile) Sterile: No Pack Size: 960 Tips (10 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p> | 960 Tips (10 Racks × 96 Tips) | S1113-1806 |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (non-sterile), Case Sterile: No Pack Size: 7,680 Tips (80 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p> | 7,680 Tips (80 Racks × 96 Tips) | S1113-1806-C |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p> | 960 Tips (10 Racks × 96 Tips) | S1113-1816 |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Racks (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Racks × 96 Tips)</p> <p> Create PDF Data Sheet</p> | 7,680 Tips (80 Racks × 96 Tips) | S1113-1816-C |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)</p> <p> Create PDF Data Sheet</p> | 960 Tips (10 Refills × 96 Tips) | S1113-1706 |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)</p> <p> Create PDF Data Sheet</p> | 9,600 Tips (100 Refills × 96 Tips) | S1113-1706-C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk







| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|--|--------------|
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (sterile) Sterile: Yes Pack Size: 960 Tips (10 Refills × 96 Tips)</p> <p> Create PDF Data Sheet</p> | 960 Tips (10 Refills × 96 Tips) | S1113-1716 |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Refills (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Refills × 96 Tips)</p> <p> Create PDF Data Sheet</p> | 7,680 Tips (80 Refills × 96 Tips) | S1113-1716-C |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack rack (non-sterile) Sterile: No Pack Size: 960 Tips (1 Stack rack × 960 Tips)</p> <p> Create PDF Data Sheet</p> | 960 Tips (1 Stack rack × 960 Tips) | S1113-1206 |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack racks (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)</p> <p> Create PDF Data Sheet</p> | 9,600 Tips (10 Stack racks × 960 Tips) | S1113-1206-C |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack rack (sterile) Sterile: Yes Pack Size: 960 Tips (1 Stack rack × 960 Tips)</p> <p> Create PDF Data Sheet</p> | 960 Tips (1 Stack rack × 960 Tips) | S1113-1216 |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Stack racks (sterile), Case Sterile: Yes Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)</p> <p> Create PDF Data Sheet</p> | 9,600 Tips (10 Stack racks × 960 Tips) | S1113-1216-C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|-----------------------------------|--------------|
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Bag (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bag × 1000 Tips)</p> <p> Create PDF Data Sheet</p> | 1,000 Tips (1 Bag × 1000 Tips) | S1113-1006 |
|  <p>200 µl UltraPoint® graduated TipOne® Tip, Yellow, Bags (non-sterile), Case Sterile: No Pack Size: 25,000 Tips (25 Bags × 1000 Tips)</p> <p> Create PDF Data Sheet</p> | 25,000 Tips (25 Bags × 1000 Tips) | S1113-1006-C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk