



10/20 µl XL graduated TipOne® Filter Tip (sterile)



Product Highlights

- Graduated at 1.0 µl, 5.0 µl, 10 µl and 20 µl
- Highly effective HDPE filter that contains no additives to interfere with samples
- Universal tip fit
- High-quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and Certified RNase, DNase, DNA and Pyrogen free
- All components, excluding tips and shrink wrap, are recyclable

Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume range	0.1 - 20 µl
Volume	20 µl
Tip Type	XL
Material	Polypropylene
Graduations	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 10/20 µl XL graduated TipOne® Filter Tip (sterile)

Experience effortless precision with TipOne® [Filter Tips](#)

Elevate your daily lab work with TipOne® [filter tips](#)! Enjoy precise, reliable pipetting every time—no more wasted samples or frustrating fit issues. TipOne® [filter tips](#) are compatible with all popular [pipettes](#) for accurate results you can trust, every single day.

Why TipOne® makes labs better:

Universal tip fit: Simplify your bench and inventory. The tip collar of TipOne® tips is engineered to fit all popular [pipettes](#) without any loss of accuracy or precision.

Made for precision: Expert design and the use of high-quality virgin materials means our tips provide consistent accuracy to help you get reliable results and reduce costly rework.

Excellent filter: The TipOne® filter protects against aerosol cross-contamination without trapping your valuable samples. And, the filters contain no additives to interfere with samples either!

Pure quality: Sterile TipOne® tips are validated according to EN ISO 11137-2:2007. All TipOne® tips are certified RNase, DNase, DNA and pyrogen free.

Flexible Packaging:

Racked: Sturdy, reusable racks stay stable on your bench. Quickly reload with our sterile or non-sterile TipOne® refills for reduced plastic waste and cost.

Refill: Load sterile filter tips into racks in an instant, and cut down on plastic waste and costs. TipOne® sterile refills are fully enclosed so tip sterility is never compromised.








Join our journey to sustainable pipetting!

Reduce your lab's plastic waste! TipOne® racks and refill system components use recycled materials and are fully recyclable (except shrink wrap and the tips themselves) via the TipOne® recycling service.

We'll collect your used TipOne® plastic, recycle it and use the recycled plastic to make more TipOne® products. No more guessing if or how your recycling gets reused—TipOne® closes the loop for you!



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 10/20 µl XL graduated TipOne® Filter Tip, Refills (sterile) Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1120-3710
 10/20 µl XL graduated TipOne® Filter Tip, Refills (sterile), Case Pack Size: 7,680 Tips (80 Refills × 96 Tips)	7,680 Tips (80 Refills × 96 Tips)	S1120-3710-C
 10/20 µl XL graduated TipOne® Filter Tip, Racks (sterile) Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1120-3810
 10/20 µl XL graduated TipOne® Filter Tip, Racks (sterile), Case Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1120-3810-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk