



Distel Dosing System



Product Highlights

- The new-style disinfectant system to replace wipes
- Effective against bacteria, yeast and viruses
- 1 L of concentrate provides 2000 x 5.0 ml doses
- Volume markers on bottles aid quick and easy dilution
- Dosing heads provide consistent dosage, every time
- Save plastic waste! Reuse the bottles again and again

Product Applications

Disinfection
Decontamination
Cleaning
Occupational health and safety

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
----------	----------------

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk





More informations about Distel Dosing System

The new-style disinfectant system to replace wipes! Take advantage of this new system where you can use any cloth/towel and dose it with the amount you need. Distel DS is a convenient and easy to use surface disinfectant for the cleaning and disinfection of hard surfaces and non-invasive medical devices. Simply dilute 1 litre of concentrate with tap water to make 10 litres of disinfectant that is effective against bacteria, yeast and viruses in just 3 minutes! The clever dosing head on the bottles ensures a consistent 5.0 ml dose of disinfectant every time. A starter pack is available that contains all you need to get going (1L disinfectant concentrate and empty reusable dosing bottles); you can then purchase the 1L concentrate on its own as and when needed. It's easy, green and results in less plastic waste. Eliminate single-use plastic wipes now with Distel DS - just dilute, dose and wipe using a towel or cloth of your choice.



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1L Distel DS Disinfectant Concentrate Packaging Unit: Bottle  Create PDF Data Sheet	1 Piece (1 Bottle × 1 Piece)	CCH050201

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk