



Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 250 ml



Product Highlights

- 0.22 µm PES membrane
- High speed filtration
- Funnels compatible with all common bottles (GL 45 compatible)
- Temperature range: -20°C to +50 C°
- Sterile SAL10⁻⁶ (ISO11137) / 3-year shelf life
- Certified RNase, DNase free and non-pyrogenic
- Wide bottle base for stability
- Wide funnel top for fast and efficient use
- Hose connector designed to fit multiple hose diameters

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

| | |
|----------------------|---------------------------|
| Art. No. | CC8221-4226 |
| Pack Size | 12 Pcs. (1 Box × 12 Pcs.) |
| Sterile | Yes |
| Volume | 250 ml |
| Free of endotoxins | Yes |
| Free of heavy metals | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Non-cytotoxic | Yes |
| Pore size | 0.22 µm |
| System | Full Assembly |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk















More informations about Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 250 ml

OPEN WIDE!

Bottle Top Filtration units with a wide base for stability and an easily accessible, wide funnel top for efficient use. The clever hose connector is designed to fit multiple hose diameters. The media bottle caps have large knurls for easy loosening or tightening. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number and expiry date for easy traceability.







All Variations

| PRODUCT NAME | | PACKAGING SIZE | ART. NO. |
|---|--|---------------------------|-------------|
|  | Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 150 ml Pore size: 0.22 µm Volume: 150 ml | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-3226 |
| |  Create PDF Data Sheet | | |
|  | Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 250 ml Pore size: 0.22 µm Volume: 250 ml | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-4226 |
| |  Create PDF Data Sheet | | |
|  | Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 500 ml Pore size: 0.22 µm Volume: 500 ml | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-5226 |
| |  Create PDF Data Sheet | | |
|  | Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 1000 ml Pore size: 0.22 µm Volume: 1000 ml | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-6226 |
| |  Create PDF Data Sheet | | |
|  | Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 150 ml Pore size: 0.45 µm Volume: 150 ml | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-3456 |
| |  Create PDF Data Sheet | | |
|  | Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 250 ml Pore size: 0.45 µm Volume: 250 ml | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-4456 |
| |  Create PDF Data Sheet | | |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



| PRODUCT NAME | | PACKAGING SIZE | ART. NO. |
|---|--|---------------------------|-------------|
|  | Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 500 ml Pore size: 0.45 µm Volume: 500 ml | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-5456 |
|  Create PDF Data Sheet | | | |
|  | Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 1000 ml Pore size: 0.45 µm Volume: 1000 ml | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-6456 |
|  Create PDF Data Sheet | | | |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk