



Bottle Top Filtration Units (sterile)



Product Highlights

- Pore sizes: 0.22 and 0.45 μm
- Compatible with standard GL 45 thread of standard bottles
- Temperature range: -20°C to $+50^{\circ}\text{C}$
- Sterile SAL10⁻⁶ (ISO11137)
- Certified DNase/RNase free and non-pyrogenic
- Certified non-cytotoxic and heavy metal free
- Reduced plastic vacuum packaging
- Crystal clear polystyrene funnel top and bottle base

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|----------------------|---------------------------|
| Art. No. | See variations |
| Sterile | Yes |
| Pack Size | 12 Pcs. (1 Box × 12 Pcs.) |
| Free of endotoxins | Yes |
| Free of heavy metals | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Non-cytotoxic | Yes |
| System | Full Assembly |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about **Bottle Top Filtration Units (sterile)**

OPEN WIDE!













[Bottle Top Filtration](#) units with a wide base for stability and a wide, easily accessible bottle top for efficient use, and the reservoir bottle caps have large knurls for easy loosening or tightening. The clever hose connector is designed to fit multiple hose diameters. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number for easy traceability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk







All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------|-------------|
|  <p>Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 150 ml Pore size: 0.22 µm Volume: 150 ml</p> <p> Create PDF Data Sheet</p> | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-3226 |
|  <p>Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 250 ml Pore size: 0.22 µm Volume: 250 ml</p> <p> Create PDF Data Sheet</p> | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-4226 |
|  <p>Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 500 ml Pore size: 0.22 µm Volume: 500 ml</p> <p> Create PDF Data Sheet</p> | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-5226 |
|  <p>Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 1000 ml Pore size: 0.22 µm Volume: 1000 ml</p> <p> Create PDF Data Sheet</p> | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-6226 |
|  <p>Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 150 ml Pore size: 0.45 µm Volume: 150 ml</p> <p> Create PDF Data Sheet</p> | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-3456 |
|  <p>Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 250 ml Pore size: 0.45 µm Volume: 250 ml</p> <p> Create PDF Data Sheet</p> | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-4456 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------|-------------|
|  Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 500 ml Pore size: 0.45 µm Volume: 500 ml  Create PDF Data Sheet | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-5456 |
|  Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 1000 ml Pore size: 0.45 µm Volume: 1000 ml  Create PDF Data Sheet | 12 Pcs. (1 Box × 12 Pcs.) | CC8221-6456 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk