

## 200 µl UltraPoint® graduated TipOne® Filter, Natural, Refills (sterile)



### Product Highlights

- Graduated at 50 µl, 100 µl and 200 µl
- UltraPoint® thin-wall orifice helps prevent sample hang up
- Highly effective HDPE filter that contains no additive to interfere with your samples
- Universal tip fit
- High-quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and Certified RNase, DNase, DNA and Pyrogen free
- Refills dramatically reduce plastic waste compared to pre-racked tips
- Refill top and bottom shells made from 100 % recycled plastic
- All components, excluding tips and shrink wrap, are recyclable

## Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|                           |                                                                                                                                                                                                                  |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Art. No.                  | S1123-8710                                                                                                                                                                                                       |
| Pack Size                 | 960 Tips (10 Refills × 96 Tips)                                                                                                                                                                                  |
| Color Code Wafer          | Yellow                                                                                                                                                                                                           |
| Sterile                   | Yes                                                                                                                                                                                                              |
| Autoclavable              | No                                                                                                                                                                                                               |
| Filter tip                | Filter                                                                                                                                                                                                           |
| Volume range              | 1.0 - 200 µl                                                                                                                                                                                                     |
| Volume                    | 200 µl                                                                                                                                                                                                           |
| Tip Type                  | UltraPoint®                                                                                                                                                                                                      |
| Material                  | Polypropylene                                                                                                                                                                                                    |
| Graduations               | Yes                                                                                                                                                                                                              |
| Color                     | Natural                                                                                                                                                                                                          |
| Sustainability attributes | Refills dramatically reduce plastic waste compared to pre-racked tips.<br>Refill top and bottom shells made from 100 % recycled plastic.<br>All components, excluding tips and shrink wrap wrap, are recyclable. |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



For use with

**Abimed:** BioPette-E, Discovery Comfort, LABMATE Optima

**Biohit:** eLine, ePet, mLine, Picus, Proline

**BRAND:** Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S

**Capp:** CappAero Comfort, CappAero Multi, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel

**Eppendorf:** Reference, Reference 2, Research, Research plus, Research pro, Xplorer

**Gilson:** Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

**Hamilton:** SoftGrip Pipette, SofTouch

**HTL:** Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

**Nichiryo:** 00-NP7, 50F, 50V, NLE, NPLO, NPX

**Rainin:** Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

**Sartorius:** eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta


**Socorex:** Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital

**STARLAB:** ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E

**Thermo Fisher Scientific Inc.:** Finn Focus Micro, Finn Focus Univ., Finnpipette Digital, Finnpipette Novus Turquoise

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 200 µl UltraPoint® graduated TipOne® Filter, Natural, Refills (sterile)

### Nothing goes through!

Protect your samples from unwanted contamination with TipOne® filter tips. Our high-quality filters prevent aerosol cross-contamination without trapping your valuable samples. Plus, they contain no additives that could interfere with your samples.

You can't always see fine aerosols entering your pipette, but they can contaminate other samples later. TipOne® filter tips block over 99% of aerosols\*, making them one of the best filters available.

Sterile TipOne® tips are validated according to EN ISO 11137-2:2007.

The filter tip refill, is an economical, plastic-saving alternative to racked filter tips.


### Reduce single-use plastic!

Compared to pre-racked tips, using TipOne® sterile filter refills dramatically reduces the amount of single-use plastic going through your laboratory. The sterile refill is quick and easy to use, just unwrap and place the entire refill into an empty TipOne® rack. The fully enclosed refill means sterility is never compromised. Tip info is printed on every wafer, for full traceability when racked.

### Join our journey to sustainable pipetting

TipOne® tip racks are built to last and can be reused again and again with these cost-effective refills. The refill's top and bottom shells are made from 100% recycled plastic. All components of TipOne® sterile filter tips refills are recyclable (excluding tips and shrink wrap). Sign up for our TipOne® recycling service!

\*See TipOne® brochure in downloads section for details.



## All Variations

| PRODUCT NAME                                                                                                                                                                                                                                                                                                                                                         | PACKAGING SIZE                    | ART. NO.     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|--------------|
| <br><a href="#">200 µl UltraPoint® graduated TipOne® Filter, Natural, Refills (sterile)</a><br><b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)<br> <a href="#">Create PDF Data Sheet</a>         | 960 Tips (10 Refills × 96 Tips)   | S1123-8710   |
| <br><a href="#">200 µl UltraPoint® graduated TipOne® Filter, Natural, Refills (sterile), Case</a><br><b>Pack Size:</b> 7,680 Tips (80 Refills × 96 Tips)<br> <a href="#">Create PDF Data Sheet</a> | 7,680 Tips (80 Refills × 96 Tips) | S1123-8710-C |
| <br><a href="#">200 µl UltraPoint® graduated TipOne® Filter, Natural, Rack (sterile)</a><br><b>Pack Size:</b> 960 Tips (10 Rack × 96 Tips)<br> <a href="#">Create PDF Data Sheet</a>              | 960 Tips (10 Rack × 96 Tips)      | S1123-8810   |
| <br><a href="#">200 µl UltraPoint® graduated TipOne® Filter, Natural, Racks (sterile), Case</a><br><b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)<br> <a href="#">Create PDF Data Sheet</a> | 7,680 Tips (80 Racks × 96 Tips)   | S1123-8810-C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)